COLUMN USER GUIDE

for Agilent InfinityLab Poroshell Chiral Columns

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This booklet provides general information for Agilent InfinityLab Poroshell 120 Chiral columns.

For additional detailed information about your specific phase or family, see:

www.agilent.com/chem/poroshell-120-chiral

Getting Started

A QC Column Performance Report, including a test chromatogram, is enclosed with every Agilent column. The QC test system has been modified from a standard system to minimize system dead volume, so it may vary from the system used in your lab. This allows a better evaluation of the column and assures a more consistent product. A properly configured LC system will generate similar results to the chromatogram on your QC Performance Report.

Modern columns are robust and are designed to operate for long periods under normal chromatographic conditions. You can maximize column performance by running it within specifications. Always review the specifications before putting in place a final method.

Using Your Column

Installation

The direction of flow is marked on the column.

For removable, zero-dead volume column connections, Agilent recommends the use of the InfinityLab Quick Connect and Quick Turn family. Choices are:

Maximum System Pressure	Recommended Fitting	Part Numbers
Up to 400 bar	InfinityLab Quick Turn Fitting (finger-tight)	Fitting: 5067-5966
Up to 800 bar	InfinityLab Quick Turn Fitting (with Mounting tool)	Fitting: 5067-5966 Mounting Tool: 5043-0915
Up to 1300 bar	InfinityLab Quick Connect Fitting	Fitting: 5067-5965

For more information and part numbers, please see to the Agilent InfinityLab Fitting Brochure (5991-5164EN).

InfinityLab Quick Connect assembly, p/n 5067-5961

InfinityLab Quick Turn fitting, p/n 5067-5966

Learn more at www.agilent.com/chem/infinitylabfittings

Column conditioning

Every column is tested before shipment and shipped in solvent, typically methanol. Care should be taken to not exceed pressure limits when transitioning mobile phase compositions. Flushing with 10 to 20 column volumes should help in transitioning to your mobile phase. The column should be properly equilibrated prior to use. This will ensure reproducibility and help prevent retention time drifting.

Important safety considerations

- All points of connection in liquid chromatographic systems are potential sources of leaks. Users should be aware of the toxicity or flammability of their mobile phases.
- Because of the small particle size, dry column packings are respirable. Columns should only be opened in a well-ventilated area.
- Please adhere to operating pressure limits noted for each column (see chart). Exceeding these limits will compromise chromatographic performance and could be unsafe.

Other operating tips

- While generally not harmful to the column, reverse flow should be avoided.
- Always use high purity reagents and chromatography grade solvent to prepare your mobile phase. Degas and filter all mobile phase prior to use.
- Disassembling a column will degrade column performance.
- New columns contain organic solvents. See your Column QC Performance Report for the solvent composition in your column. Initially, care should be taken not to pass any mobile phase through the column that may cause a precipitate to form.
- Agilent chiral columns are compatible with water and all common organic solvents.
- Columns should not be maintained at elevated pH or elevated temperature when not in use.
- Avoid use of this column outside of recommended pressure, temperature, and pH ranges for column phase (see next page). Expect reduced lifetime when operating outside the recommended ranges.

Column Operating Parameters: pH and Temperature

Column	Chiral selector (covalently bonded onto Agilent Poroshell particles)	Max operating temperature	Recommended pH range
Agilent InfinityLab Poroshell Chiral-T	Teicoplanin	45	2.5-7
Agilent InfinityLab Poroshell Chiral-V	Vancomycin	45	2.5-7
Agilent InfinityLab Poroshell Chiral-CD	Hydroxypropylated- β - cyclodextrin	45	3-7
Agilent InfinityLab Poroshell Chiral-CF	Derivatized cyclofructan (CF6)	45	3-7

Note: All silica-based packings have some solubility in pH >6 aqueous mobile phases. When using silica-based columns at pH >6, best column lifetime is obtained at lower temperatures (40 °C max) using low buffer concentrations in the range of 0.01 to 0.02 M. Operating at extreme ends of pH and temperature ranges will have a significant impact on column lifetime.

Maximum Operating Pressures

Column	Particle size	Max pressure
Agilent InfinityLab Poroshell Chiral-T	2.7 μ m	400 bar (6,000 psi)
Agilent InfinityLab Poroshell Chiral-V	2.7 μ m	400 bar (6,000 psi)
Agilent InfinityLab Poroshell Chiral-CD	2.7 μ m	400 bar (6,000 psi)
Agilent InfinityLab Poroshell Chiral-CF	2.7 μ m	400 bar (6,000 psi)

Mobile Phase Selection and Operating Temperatures

Agilent Chiral columns can be used in any type of organic solvents. However, when switching to or from a normal phase condition to other mobile phases, please use ethanol or IPA as a transition solvent. Allow at least 10 column volumes of conditioning time before injection. **Care should be taken to not exceed the pressure limit of the column when transitioning mobile phases.**

Recommended Screening Mobile Phases

Column	Screening mobile phases (isocratic)
Agilent InfinityLab Poroshell Chiral-T and Chiral-V	<p>Polar ionic: 100/0.2 wt% MeOH/ammonium formate</p> <p>Polar organic: 60/40/0.3/0.2 ACN/MeOH/acetic acid/TEA</p> <p>Reversed-phase: 30/70 MeOH/20 mM ammonium formate (pH 4)</p> <p>Normal-phase: 20/80/0.3/0.2 EtOH/heptane/acetic acid/TEA</p>
Agilent InfinityLab Poroshell Chiral-CD	<p>Reversed-phase: 30/70 ACN/20 mM Ammonium formate buffer (pH = 4)</p> <p>Polar organic: 60/40/0.3/0.2 ACN/MeOH/acetic acid/TEA</p>
Agilent InfinityLab Poroshell Chiral-CF	<p>Polar organic: 60/40/0.3/0.2 ACN/MeOH/acetic acid/TEA</p> <p>Normal-phase: 20/80/0.3/0.2 EtOH/heptane/acetic acid/TEA</p>

Recommended Method Optimizations

Column	Recommended method optimization
Agilent InfinityLab Poroshell Chiral-T and Chiral-V	<p>Polar ionic mode: Change concentration of salts between 0.01 wt% and 1 wt%</p> <p>Polar organic: Change ratio of ACN to MeOH</p> <p>Reversed-phase: Change organic/buffer ratio Change different organic solvents (ACN or MeOH) Change pH for ionizable compounds (typically pH 3-5 for bases and pH 5-7 for acids)</p> <p>Normal-phase: Change ratio of alcohol to alkanes, change type of alcohols</p>
Agilent InfinityLab Poroshell Chiral-CD	<p>Reversed-phase: Change organic/buffer ratio Change pH for ionizable compounds</p> <p>Polar organic: Change ratio of ACN to MeOH</p>
Agilent InfinityLab Poroshell Chiral-CF	<p>Polar organic: Change ratio of ACN to MeOH</p> <p>Normal-phase: Change ratio of alcohol to alkanes, and/or change type of alcohol</p>

Column Care

Column regeneration

Column	Regeneration procedure
Agilent InfinityLab Poroshell Chiral-T and Chiral-V	Flush the column with 50/50, ACN/50 mM ammonium acetate at lower flow rates for at least 2 hours. Then flush with pure methanol, and equilibrate with mobile phase.
Agilent InfinityLab Poroshell Chiral-CD and Chiral-CF	Flush the column with 50/50, ACN/MeOH at lower flow rates for at least 2 hours. Then equilibrate with mobile phase.

Storage Recommendations

Long-term storage of silica-based, bonded phase columns should be in a pure organic solvent. If the column has previously been used with a buffered mobile phase, the buffer should first be removed by purging the column with 20 to 30 column volumes of a 50:50 mixture of methanol or acetonitrile and water, followed by 20 to 30 column volumes of the pure solvent. Before storing, end fittings should be tightly capped with end plugs to prevent packing from drying out. Columns may be safely stored for short periods in most mobile phases. To protect equipment, is it desirable to remove salts from the instrument and column by purging the column with the same mobile phase without the buffer (for example, using 30/70 MeOH/H₂O to remove a 30/70 MeOH/20 mM ammonium formate buffered mobile phase). Re-equilibration is rapid with the original mobile phase when using this approach and any danger of corrosion from the salts is eliminated. Refer to recommended storage solvents in the table below.

Column Storage

Column	Regeneration procedure
InfinityLab Poroshell Chiral-T and V	Pure MeOH or ethanol is recommended for long term storage.
InfinityLab Poroshell Chiral-CD and CF	Pure ethanol is recommended for long term storage.

Use of Chiral columns in Supercritical Fluid Chromatography (SFC)

Agilent Poroshell 120 Chiral columns are suitable for SFC. Column pressure and pH limits remain valid in SFC. Please keep in mind that pH in SFC is not directly comparable with pH values of aqueous solutions. Therefore, the apparent pH in CO₂/organic modifier might differ.

Column equilibration

Before use, flush the column with 10–30 volumes of methanol or isopropyl alcohol to completely flush out shipping solvents. Then, equilibrate the column with approximately 10–20 column volumes of mobile phase at starting conditions.

Column storage

Make sure that you completely remove additives such as acids, bases, or buffers to prevent damage to the column. Flush the column with 10–20 column volumes of ethanol to remove all carbon dioxide. Do not store columns in a mixture of CO₂ and MeOH, as CO₂ will evaporate and dry out the stationary phase. Properly seal with end plugs to prevent solvent evaporation.

Tips for Getting the Best Chromatographic Results

- Optimize your instrument by minimizing tubing lengths between components, to reduce extra column volume and band broadening. Use 0.12 mm id red tubing for Fast LC/high efficiency columns. Learn about capillary options at [agilent.com/chem/lccapillaries](https://www.agilent.com/chem/lccapillaries)
- Ensure the data collection rate is optimized for your column. Use a higher collection rate for Fast LC columns (InfinityLab Poroshell 120, RRHT, RRHD).
- Use sample filtration or other sample prep as appropriate for your sample. Learn more at [agilent.com/chem/sampleprep](https://www.agilent.com/chem/sampleprep)
- Use Agilent certified lamps in your LC instruments for best performance.

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
Further Information

For help choosing the right column, request the
LC Column Selection Guide or the
BioHPLC Column Selection Guide
agilent.com/chem/getguides

The LC Column and Sample Prep Navigator is a web-based tool that can help you find the right column for your method, based on your input. Consult agilent.com/chem/navigator

Agilent Ordering Information

For more information on our products and services,
visit our web site at **agilent.com**

For technical support and local information,
visit **agilent.com/chem/columnsupport**

To place an order,
visit **agilent.com/chem/wheretobuy**

Agilent offers a complete line of sample preparation products to
support LC and LC/MS applications.

The Agilent Bond Elut SPE and Captiva Filtration Sample Prep
family of products offer the widest range of solutions for every
level of sample cleanliness to help you increase throughput and
enhance the quality of your data.

Learn more at **agilent.com/chem/sampleprep**

This information is subject to change without notice.

Agilent Technologies, Inc. 2017, 2019

Printed in USA, October 8, 2019

820120-014