

Syringes | LC autosampler

Manufactured since 2004, the comprehensive range of L-MARK® Gold and L-MARK syringes for LC systems are designed and tested to meet critical instrument parameters.

Features and benefits

- Plunger guide maintains optimal plunger alignment for extended plunger lifetime.
- Syringes available with polyethylene, PTFE and black PTFE tipped plungers.
- Titanium plungers for longer lifetime.
- Needle setting in the syringe barrel is protected by a PTFE barrier for reduced exposure of the adhesive to samples.
- The L-MARK Gold syringes include a chemically inert needle for extremely low carryover.

Product specifications

Available for CTC Analytics and Thermo Scientific LC autosamplers.

Specifications	
Fixed needle syringes	
Accuracy and reproducibility	±1% (dispensed volume)
Temperature range	5-70°C

For more information about this product visit www.trajanscimed.com or contact techsupport@trajanscimed.com

Specifications are subject to change without notice.

Analytical syringes manufactured by Trajan Scientific and Medical are intended for analytical and laboratory use only and are not intended or approved for use with food, including the production or packaging of food, nor medical or human in-vivo use.

