


MACHEREY-NAGEL

Dip & Read Tests

Rapid Tests



Disinfectant testing

Rapid Tests for Healthcare

- Quick & easy
- Reliable results
- Higher patient safety



MACHEREY-NAGEL

www.mn-net.com



MN Dip & Read Tests for Healthcare

Healthcare professionals use different types of disinfectants and sanitizers for general disinfection and for cleaning of medical devices. Prior to the disinfection of devices it proved useful to control the potency of the disinfection baths. After disinfection all devices have to be rinsed and controlled for sanitizer residues prior to the next use. Dip & Read tests allow the quick and easy detection of sanitizer residues. Thus, MN test strips increase patient safety, while making healthcare professionals' work easier and more efficient.



Product information

Disinfectant	Product	Potency	Residual	Gradation	REF
Chlorine	QUANTOFIX® Chlorine 500	■		0 · 25 · 50 · 100 · 250 · 500 mg/L Cl ₂	91354
	QUANTOFIX® Chlorine Sensitive*		■	0 · 0.1 · 0.5 · 1 · 3 · 10 mg/L Cl ₂	91339
Chlorine dioxide	QUANTOFIX® Chlorine dioxide 3000	■		0 · 10 · 25 · 50 · 100 · 250 · 500 · 1000 · 1500 · 2000 · 2500 · 3000 mg/L ClO ₂	91357
	QUANTOFIX® Chlorine dioxide 15		■	0 · 0.25 · 0.5 · 1 · 3 · 5 · 10 · 15 mg/L ClO ₂	91355
Glutaraldehyde	QUANTOFIX® Glutaraldehyde	■		0 · 0.5 · 1.0 · 1.5 · 2.0 · 2.5 % glutaraldehyde	91343
Hydroxyacetic Acid	pH-Fix 3.6-6.1*		■	pH 3.6 · 4.1 · 4.4 · 4.7 · 5.0 · 5.3 · 5.6 · 6.1	92130
Peracetic Acid	QUANTOFIX® Peracetic acid 2000*	■		0 · 500 · 1000 · 1500 · 2000 mg/L peracetic acid	91342
	QUANTOFIX® Peracetic acid 500*	■		0 · 50 · 100 · 200 · 300 · 400 · 500 mg/L peracetic acid	91341
	QUANTOFIX® Peracetic acid 50*		■	0 · 5 · 10 · 20 · 30 · 50 mg/L peracetic acid	91340
Peroxide	QUANTOFIX® Peroxide 1000*	■		0 · 50 · 150 · 300 · 500 · 800 · 1000 mg/L H ₂ O ₂	91333
	QUANTOFIX® Peroxide 100*		■	0 · 1 · 3 · 10 · 30 · 100 mg/L H ₂ O ₂	91312
	QUANTOFIX® Peroxid 25*		■	0 · 0.5 · 2 · 5 · 10 · 25 mg/L H ₂ O ₂	91319
QUAT	QUANTOFIX® QUAT	■		0 · 10 · 25 · 50 · 100 · 250 · 500 · 1000 mg/L benzalkoniumchloride	91337

* suitable for evaluation with QUANTOFIX® Relax

QUANTOFIX® Relax

The instrument is for the automatic reading of QUANTOFIX® test strips. Just place the strip on the strip adaptor and the Relax will give a precise and reliable result within seconds. The results are printed, saved internally and can be transmitted to a PC. It is ideal for QC managers and decision makers in production.

QUANTOFIX® Relax – for easy, reliable QC procedures

- highest precision
- objective results
- easy result management

REF 91346



www.mn-net.com

MACHEREY-NAGEL



MACHEREY-NAGEL GmbH & Co. KG
Valencienner Str. 11
52355 Düren · Germany

DE Tel.: +49 24 21 969-0 info@mn-net.com
CH Tel.: +41 62 388 55 00 sales-ch@mn-net.com
FR Tel.: +33 388 68 22 68 sales-fr@mn-net.com
US Tel.: +1 888 321 62 24 sales-us@mn-net.com

