

# Safety Data Sheet

according to Regulations REACH 1907/2006/EC

REF: 985680

NANOCOLOR tot.Phosphate 15, robot

Page: 1/14

Printing date: 15.05.2024

Date of issue: 26.09.2022

Version: 2.2.3.10

## SECTION 1: Identification of the substance/mixture and of the company

### 1.1 Product identifier

REF 985680  
Product name NANOCOLOR tot.Phosphate 15, robot

REACH Registration number(s): see SECTION 3.1/3.2 or  
A registration number for the substance(s) does not exist because the annual tonnage does not require registration or the substance or its use is excluded from registration.

1 x 2 g total Phosphate (B1)  
1 x 5 mL total Phosphate (B2)  
20 x 35 mg NANOFIX total Phosphate 15 (R0) UFI: 9M7U-D3PX-H206-F86E  
1 x 110 mL total Phosphate 15 Reagent A  
1 x 6 mL total Phosphate (C) UFI: 7F7U-D3A4-W206-4K19

### 1.2 Relevant identified uses of the substance or mixture and uses advised against

#### Relevant identified uses

Product for analytical use.

Exposure Scenario Classification according REACH, RIP 3.2 Codes: SU 0-2, PC 21, PROC 15, AC 0  
The exposure scenario is integrated into sections 1-16.

#### Uses advised against

not described

### 1.3 Details of the supplier of the safety data sheet

#### Manufactured by:

MACHEREY-NAGEL GmbH & Co. KG  
Valenciener Str. 11, 52355 Düren, Germany  
Phone: +49 2421 969 0

E-mail: [sds@mn-net.com](mailto:sds@mn-net.com) ([msds@mn-net.com](mailto:msds@mn-net.com))

### 1.4 Emergency telephone number

Outside Germany (DE): Call your regional Poisons Information Service or call local Life Saving Service.

DE: Gemeinsames Gif tinformat ionszentrum (GGIZ)

99089 Erfurt tel. +49 361 730 730, <<https://www.ggiz-erfurt.de>>

You find our current versions of SDS in Internet:

<<http://www.mn-net.com/SDS>>

## SECTION 2: Hazard identification

### 2.0 Classification of the complete product according to Regulation (EC) 1272/2008



GHS03 GHS07 GHS08

Signal word DANGER

Hazard identification	Hazard classes/categories
H272	Ox. Liq. 2
H302	Acute Tox. 4 oral
H315	Skin Irrit. 2
H317	Skin Sens. 1
H319	Eye Irrit. 2
H334	Resp. Sens. 1
H335	resp. irrit. STOT SE 3

### 2.1 Classification of the substance or mixture according to Regulation (EC) 1272/2008

6 mL total Phosphate (C)



MACHEREY-NAGEL GmbH & Co. KG  
Valenciener Str. 11  
52355 Düren · Germany  
[www.mn-net.com](http://www.mn-net.com)

DE Tel.: +49 24 21 969-0 [info@mn-net.com](mailto:info@mn-net.com)  
CH Tel.: +41 62 388 55 00 [sales-ch@mn-net.com](mailto:sales-ch@mn-net.com)  
FR Tel.: +33 388 68 22 68 [sales-fr@mn-net.com](mailto:sales-fr@mn-net.com)  
US Tel.: +1 888 321 62 24 [sales-us@mn-net.com](mailto:sales-us@mn-net.com)

# Safety Data Sheet

according to Regulations REACH 1907/2006/EC

REF: 985680

NANOCOLOR tot.Phosphate 15, robot

Page: 2/14

Printing date: 15.05.2024

Date of issue: 26.09.2022

Version: 2.2.3.10



GHS07

Signal word                      WARNING

**Hazard identification**                      **Hazard classes/categories**

H315                                      Skin Irrit. 2  
 H319                                      Eye Irrit. 2

**2 g total Phosphate (B1)**

Signal word                                      Do not need labelling as hazardous  
 -

No hazard class

**110 mL total Phosphate 15 Reagent A**

Signal word                                      Do not need labelling as hazardous  
 -

No hazard class

**5 mL total Phosphate (B2)**

Signal word                                      Do not need labelling as hazardous  
 -

No hazard class

**35 mg NANOFIX total Phosphate 15 (R0)**



GHS03



GHS07



GHS08

Signal word                                      DANGER

**Hazard identification**                      **Hazard classes/categories**

H272                                      Ox. Liq. 2  
 H302                                      Acute Tox. 4 oral  
 H315                                      Skin Irrit. 2  
 H317                                      Skin Sens. 1  
 H319                                      Eye Irrit. 2  
 H334                                      Resp. Sens. 1  
 H335                                      resp. irrit. STOT SE 3

List of H phrases: see section 16.2

**2.2 Label elements according regulation (EC) 1272/2008**

According **CLP directive** inner packages must be only labelled with GHS symbol(s) and product identifier(s) (EU 1272/2008 Annex I - 1.5.1.2). Inner packages up to 10 mL need max. 2 symbols (Annex I - 1.5.2.4.1 / 2). Harmful chemicals/mixtures with signal word: **WARNING** must not be labelled with H and P phrases **until 125 mL** (EU 1272/2008 Annex I - 1.5.2). This labelling exemption is NOT valid for sensiblizing substances. Oxidizing mixtures with signal word: **DANGER** and **H272** must not be labelled with H and P phrases **until 125 mL** .

**6 mL total Phosphate (C)**



MACHEREY-NAGEL GmbH & Co. KG  
 Valencienner Str. 11  
 52355 Düren · Germany  
[www.mn-net.com](http://www.mn-net.com)

DE Tel.: +49 24 21 969-0 info@mn-net.com  
 CH Tel.: +41 62 388 55 00 sales-ch@mn-net.com  
 FR Tel.: +33 388 68 22 68 sales-fr@mn-net.com  
 US Tel.: +1 888 321 62 24 sales-us@mn-net.com

# Safety Data Sheet

according to Regulations REACH 1907/2006/EC

REF: 985680

NANOCOLOR tot.Phosphate 15, robot

Page: 3/14

Printing date: 15.05.2024

Date of issue: 26.09.2022

Version: 2.2.3.10



GHS07

Signal word: WARNING

**2 g total Phosphate (B1)**

Do not need labelling as hazardous

Signal word: -

**110 mL total Phosphate 15 Reagent A**

Do not need labelling as hazardous

Signal word: -

**5 mL total Phosphate (B2)**

Do not need labelling as hazardous

Signal word: -

**35 mg NANOFIX total Phosphate 15 (R0)**



GHS03



GHS07



GHS08

Signal word: DANGER

H317, H334

May cause an allergic skin reaction. May cause allergy or asthma symptoms or breathing difficulties if inhaled.

P261sh, P280sh, P284, P302+352, P333+313, P342+311, P362+364, P501

Avoid breathing dust/vapours. Wear protective gloves/eye protection. [In case of inadequate ventilation] wear respiratory protection. IF ON SKIN: Wash with plenty of water. If skin irritation or rash occurs: Get medical advice/attention. If experiencing respiratory symptoms: Call a POISON CENTER/doctor. Take off contaminated clothing and wash it before reuse. Dispose of contents/container to regulated waste treatment.

### Label elements of the complete product



GHS03



GHS07



GHS08

Signal word: DANGER

H317, H334

May cause an allergic skin reaction. May cause allergy or asthma symptoms or breathing difficulties if inhaled.

P261sh, P280sh, P284, P302+352, P333+313, P342+311, P362+364, P501

Avoid breathing dust/vapours. Wear protective gloves/eye protection. [In case of inadequate ventilation] wear respiratory protection. IF ON SKIN: Wash with plenty of water. If skin irritation or rash occurs: Get medical advice/attention. If experiencing respiratory symptoms: Call a POISON CENTER/doctor. Take off contaminated clothing and wash it before reuse. Dispose of contents/container to regulated waste treatment.

## 2.3 Other hazards

### Possible hazards from physicochemical properties

In the case of pH values are less than 5 or higher than 9 then it is irritant.

### Information pertaining to particular risks to human and possible symptoms

Cause after oral intake, inhalation of vapours/dust, skin contact, impairments of health when ingested in small quantities. May cause sensitization by skin contact, also in repeated contact of small amounts. May cause allergy or asthma symptoms or breathing difficulties if inhaled.

### Information pertaining to particular risks to the environment

PBT: not applicable

vPvB: not applicable

### Possible endocrine disrupting effects



MACHEREY-NAGEL GmbH & Co. KG  
Valenciener Str. 11  
52355 Düren · Germany  
[www.mn-net.com](http://www.mn-net.com)

DE Tel.: +49 24 21 969-0 [info@mn-net.com](mailto:info@mn-net.com)  
CH Tel.: +41 62 388 55 00 [sales-ch@mn-net.com](mailto:sales-ch@mn-net.com)  
FR Tel.: +33 388 68 22 68 [sales-fr@mn-net.com](mailto:sales-fr@mn-net.com)  
US Tel.: +1 888 321 62 24 [sales-us@mn-net.com](mailto:sales-us@mn-net.com)

# Safety Data Sheet

according to Regulations REACH 1907/2006/EC

REF: 985680	NANOCOLOR tot.Phosphate 15, robot	Page: 4/14
Printing date: 15.05.2024	Date of issue: 26.09.2022	Version: 2.2.3.10

no data available

## SECTION 3: Composition / information on ingredients

### 3.1 Substances or 3.2 Mixtures

#### 6 mL total Phosphate (C)

Substance name: *ammonium heptamolybdate*  
 CAS No.: 12054-85-2

Substance rating: No criteria for classification or naming of chemical not required.  
 Formula:  $H_{24}Mo_7N_6O_{24}$   
 Pseudonym (de): Ammoniummolybdat  
 REACH Reg. No.: 01-2119498057-28-xxxx  
 EC No.: 234-722-4  
 Concentration: 2 - <5 % Correlation factor: x 0.58 (= %Mo)  
 The classification refers to the weight percentage of the metal (according to CLP regulation 2008/1272/EG Annex VI, 1.1.3.2 Note 1)  
 acc. CLP (GHS): The criteria for classification are not fulfilled.

Substance name: *sulfuric acid*  
 CAS No.: 7664-93-9

Substance rating: H315, Skin Irrit. 2, H319, Eye Irrit. 2  
 Formula:  $H_2SO_4 \cdot H_2O$   
 REACH Reg. No.: 01-2119458838-20-xxxx  
 EC No.: 231-639-5 Indice No.: 016-020-00-8  
 Concentration: 5 - <15 %  
 acc. CLP (GHS): H315, Skin Irrit. 2, H319, Eye Irrit. 2

#### 2 g total Phosphate (B1)

Substance name: *sodium L(+)-ascorbate*  
 CAS No.: 134-03-2

Substance rating: No criteria for classification or naming of chemical not required.  
 Formula:  $C_6H_7NaO_6$   
 Pseudonym (de): Ascorbinsäure, Na-Salz, E301  
 EC No.: 205-126-1  
 Concentration: 90 - <100 %  
 acc. CLP (GHS): The criteria for classification are not fulfilled.

#### 110 mL total Phosphate 15 Reagent A

Substance name: *sulfuric acid*  
 CAS No.: 7664-93-9

Substance rating: H315, Skin Irrit. 2, H319, Eye Irrit. 2  
 Formula:  $H_2SO_4 \cdot H_2O$   
 REACH Reg. No.: 01-2119458838-20-xxxx  
 EC No.: 231-639-5 Indice No.: 016-020-00-8  
 Concentration: 1 - <5 %  
 acc. CLP (GHS): The criteria for classification are not fulfilled.

#### 5 mL total Phosphate (B2)

Substance name: *water*  
 CAS No.: 7732-18-5

Substance rating: No criteria for classification or naming of chemical not required.  
 Formula:  $H_2O$   
 REACH Reg. No.: exempt, Annex IV  
 EC No.: 231-791-2  
 Concentration: 90 - <100 %  
 acc. CLP (GHS): The criteria for classification are not fulfilled.



# Safety Data Sheet

according to Regulations REACH 1907/2006/EC

REF: 985680

NANOCOLOR tot.Phosphate 15, robot

Page: 5/14

Printing date: 15.05.2024

Date of issue: 26.09.2022

Version: 2.2.3.10

## 35 mg NANOFIX total Phosphate 15 (R0)

Substance name: *sodium peroxodisulfate*  
CAS No.: 7775-27-1

Substance rating: H272, Ox. Sol. 2, H302, Acute Tox. 4 oral, H315, Skin Irrit. 2, H317, Skin Sens. 1, H319, Eye Irrit. 2, H334, Resp. Sens. 1, H335, resp. irrit. STOT SE 3

Formula:  $\text{Na}_2\text{O}_8\text{S}_2$

Pseudonym (de): Natriumpersulfat

REACH Reg. No.: 01-2119495975-15-xxxx

EC No.: 231-892-1

Concentration: 80 - <100 %

acc. CLP (GHS): H272, Ox. Liq. 2, H302, Acute Tox. 4 oral, H315, Skin Irrit. 2, H317, Skin Sens. 1, H319, Eye Irrit. 2, H334, Resp. Sens. 1, H335, resp. irrit. STOT SE 3

### 3.3 Remarks

When not listed, mixtures are added with water [CAS No. 7732-18-5] to 100%. List of H and P phrases: see section 16.2.

## SECTION 4: First aid measures

### 4.1 Description of first aid measures

Place insured person out of danger zone to fresh air immediately. Ensure quiet, warmth, and provide resuscitation if necessary. If necessary contact medical advice. Take to a doctor, in a raised position if there are breathing difficulties.

#### 4.1.1 After SKIN Contact

Remove contaminated clothing immediately. Rinse the affected skin or mucous membrane thoroughly for min. 15 minutes under running water. (If possible) use soap. Avoid neutralisation. Then apply a loose bandage.

#### 4.1.2 After EYE Contact

After contact with the eyes rinse thoroughly under running water with the eyelid wide open with eye washing bottle, eye douche or running water (protect intact eye).

#### 4.1.3 After INHALATION of vapours

After inhalation of foam or vapour fresh air should be inhaled. Keep airways free. Administer a Dexamethasone spray as soon as possible. Ensure quiet, warmth, and provide resuscitation if necessary. In the event of respiratory distress ensure that the patient inhales oxygen. Secure the breathing, heart and circulatory function. ---

#### 4.1.4 After ORAL Intake

After oral intake lots of water should be drunk after it has been ingested.

### 4.2 Most important symptoms and effects, both acute and delayed

May cause allergy or asthma symptoms or breathing difficulties if inhaled. Chronic effects: Repeated contact, even in small amounts, can lead to sensitization.

### 4.3 Indication of any immediate medical attention and special treatment needed

After SKIN CONTACT rinse with water for a long time. Apply glucocorticosteroides following inflammatory reactions. Inform patient respectively further measures and the possibility of long-term damages. ---

## SECTION 5: Firefighting measures

### 5.1 Extinguishing media

#### 5.1.1 Suitable extinguishing media

Fire extinguishers appropriate to the fire classification, and, if applicable, a fire blanket must be available in a prominent location in the work area. All extinguishers like FOAM, WATER SPRAY, DRY POWDER, CARBON DIOXIDE can be used.

#### 5.1.2 Unsuitable extinguishing media

no data available

### 5.2 Special hazards arising from the substance or mixture

Formation of hazardous and caustic vapour-air mixtures possible.

### 5.3 Advice for firefighters

No, for listed product. Product package burns like paper or plastic. Spray any vapours released with water. Retent fire water. Use only acid-resistant safety equipment.

For great amount - if necessary - protective breathing apparatus which is independent of the ambient air (isolated equipment), and sealed protective clothing is necessary in the event of a large-scale formation of toxic substances.



MACHEREY-NAGEL GmbH & Co. KG  
Valencienner Str. 11  
52355 Düren · Germany  
[www.mn-net.com](http://www.mn-net.com)

DE Tel.: +49 24 21 969-0 [info@mn-net.com](mailto:info@mn-net.com)  
CH Tel.: +41 62 388 55 00 [sales-ch@mn-net.com](mailto:sales-ch@mn-net.com)  
FR Tel.: +33 388 68 22 68 [sales-fr@mn-net.com](mailto:sales-fr@mn-net.com)  
US Tel.: +1 888 321 62 24 [sales-us@mn-net.com](mailto:sales-us@mn-net.com)

# Safety Data Sheet

according to Regulations REACH 1907/2006/EC

REF: 985680	NANOCOLOR tot.Phosphate 15, robot	Page: 6/14
Printing date: 15.05.2024	Date of issue: 26.09.2022	Version: 2.2.3.10

## 5.4 Additional information

Danger for environment **only in the event of a large-scale leakage** or formation of hazardous substances.

## SECTION 6: Accidental release measures

### 6.1 Personal precautions, protective equipment and emergency procedures

Do not breathe vapours. Wear suitable protective gloves (see 8.2.2). Wear eye protection. Regular staff training is necessary, indicating hazards and precautions on the basis of operating instructions. Restrictions on activity must be observed.

### 6.2 Environmental precautions

PBT: not applicable  
vPvB: not applicable

### 6.3 Methods and material for containment and cleaning up

Bind any escaping liquid with inert absorbent. And dispose in accordance to local regulations for the disposal of hazardous chemicals. Clean any contaminated equipment and floors with plenty of water. Collect small amounts of leaked liquid and flush with water into drains.

### 6.4 Reference to other sections

see information in section 5.4,7,8 and 13

## SECTION 7: Handling and storage

### 7.1 Precautions for safe handling

Handling in accordance with the test instruction, that comes with the product. Use only in well-ventilated working areas. Use a safety bottle when shaking test tubes.

### 7.2 Conditions for safe storage, including any incompatibilities

Safe storage is guaranteed in the original packaging from MACHEREY-NAGEL. Storage class (German chemical industry): see chapter 12.1

Storage class (VCI): 5.1B  
Water hazard class (DE): 1

#### 7.2.1 Requirements for stock rooms and containers

Keep original product packages tightly closed during handling and storage.

### 7.3 Specific end use(s)

Product for analytical use.

## SECTION 8: Exposure controls /personal protection

### 8.1 Control parameters

#### 35 mg NANOFIX total Phosphate 15 (R0)

Chemical: *sodium peroxodisulfate* CAS No.: 7775-27-1  
NIOSH: not listed  
[TWA] Time-weighted average to a reference period of 8 hours, [STEL] Short-term exposure limit related to a 15-minute period  
OSHA: not listed

#### 6 mL total Phosphate (C)

Chemical: *ammonium heptamolybdate* CAS No.: 12054-85-2  
TRGS 900 (DE): [Mo] 5 E mg/m<sup>3</sup>  
E/e respirable  
SUVA(CH) MAK value: [Mo] 5 e mg/m<sup>3</sup>

#### Chemical: *sulfuric acid* CAS No.: 7664-93-9

DNEL: 50 µg/m<sup>3</sup>  
DNEL = Derived No-Effect Level (for workers)  
PNEC (fresh water): 2.5 µg/L  
PNEC = Predicted No Effect Concentration  
TRGS 900 (DE): 0.1 E mg/m<sup>3</sup>  
E/e respirable  
Short-term exposure factor: 1 (I)  
skin resorptive (H), respiratory sensitizable (Sa), skin sensitizable (Sh), teratogenic (Z) not securely excluded / (Y) certainly excluded  
SUVA(CH) MAK value: 0,1 e mg/m<sup>3</sup>  
NIOSH: NTP Report on Carcinogens (RoC) List Yes (Known to be a human carcinogen); TWA 1 mg/m<sup>3</sup>  
[TWA] Time-weighted average to a reference period of 8 hours, [STEL] Short-term exposure limit related to a 15-minute period



# Safety Data Sheet

according to Regulations REACH 1907/2006/EC

REF: 985680	NANOCOLOR tot.Phosphate 15, robot	Page: 7/14
Printing date: 15.05.2024	Date of issue: 26.09.2022	Version: 2.2.3.10

OSHA: [TWA] 1 mg/m<sup>3</sup>

**110 mL total Phosphate 15 Reagent A**

Chemical: *sulfuric acid* CAS No.: 7664-93-9

DNEL: 50 µg/m<sup>3</sup>  
DNEL = Derived No-Effect Level (for workers)

PNEC (fresh water): 2.5 µg/L  
PNEC = Predicted No Effect Concentration

TRGS 900 (DE): 0.1 E mg/m<sup>3</sup>  
E/e respirable

Short-term exposure factor: 1 (I)  
skin resorptive (H), respiratory sensitizable (Sa), skin sensitizable (Sh), teratogenic (Z) not securely excluded / (Y) certainly excluded

SUVA(CH) MAK value: 0,1 e mg/m<sup>3</sup>  
 NIOSH: NTP Report on Carcinogens (RoC) List Yes (Known to be a human carcinogen); TWA 1 mg/m<sup>3</sup>  
[TWA] Time-weighted average to a reference period of 8 hours, [STEL] Short-term exposure limit related to a 15-minute period

OSHA: [TWA] 1 mg/m<sup>3</sup>

**2 g total Phosphate (B1)**

Chemical: *sodium L(+)-ascorbate* CAS No.: 134-03-2

**5 mL total Phosphate (B2)**

Chemical: *water* CAS No.: 7732-18-5

## 8.2 Exposure controls

Good ventilation and extraction system in the room, floor resistant to chemicals with floor drainage and washing facilities. The highest level of cleanliness must be maintained at the workplace.

### 8.2.1 Respiratory protection

Use for open access of these substances for example a protection filter, class A/AX. No additional recommendations.

### 8.2.2 Skin protection / Hand protection

Yes, gloves according EN 374 (permeation time >30 min - level 2), consist of PVC, natural latex, Neopren, or Nitril (f.ex. from Ansell or KCL). Use for short times chemical resistant latex gloves with code EN 374-3 level 1.

### 8.2.3 Eye / Face Protection

Yes, safety glasses according EN 166 with integrated side shields or wrap-around protection.

### 8.2.4 Skin protection

Recommended to avoid contamination with these hazards.

### 8.2.5 Personal hygiene

Eating, drinking, smoking, taking snuff and storage of food in work areas and at outdoor workplaces is prohibited. Avoid contact with the skin, eyes and clothing. Rinse any clothing on which the substance has been spilled, and soak it in water. Wash hands thoroughly with soap and water when stopping work and before eating, and then apply protective skin cream.

### 8.2.6 Thermal hazards

no data available

## 8.3 Limitation and monitoring of environmental exposure

Do not release product into environment.

# SECTION 9: Physical and chemical properties

## 9.1 Information on basic physical and chemical properties

**35 mg NANOFIX total Phosphate 15 (R0)**

a) State of aggregation:	solid (lyophilized)
b) Colour:	colourless
c) Odor:	odorless
d) Melting point:	no data available
e) Boiling point:	no data available
f) Flammability:	no data available
g) Explosive limits (lower / upper):	no data available
h) Flash point:	no data available
i) Flashing temperature:	no data available
j) Decomposition temperature:	no data available
k) pH value:	5-7
l) Kinematic viscosity:	no data available
m) Solubility in water:	0-100 %



MACHEREY-NAGEL GmbH & Co. KG  
 Valenciener Str. 11  
 52355 Düren · Germany  
[www.mn-net.com](http://www.mn-net.com)

DE Tel.: +49 24 21 969-0 info@mn-net.com  
 CH Tel.: +41 62 388 55 00 sales-ch@mn-net.com  
 FR Tel.: +33 388 68 22 68 sales-fr@mn-net.com  
 US Tel.: +1 888 321 62 24 sales-us@mn-net.com

# Safety Data Sheet

according to Regulations REACH 1907/2006/EC

REF: 985680

NANOCOLOR tot.Phosphate 15, robot

Page: 8/14

Printing date: 15.05.2024

Date of issue: 26.09.2022

Version: 2.2.3.10

- n) Dispersion coefficient ( $K_{o/w}$ ): no data available
- o) Vapour pressure (20°C): no data available
- p) Specific gravity: no data available
- q) Relative vapour density ( $air=1$ ): no data available
- r) Particle size: no data available

**6 mL total Phosphate (C)**

- a) State of aggregation: liquid
- b) Colour: slightly yellow
- c) Odor: odorless
- d) Melting point: no data available
- e) Boiling point: no data available
- f) Flammability: no data available
- g) Explosive limits (lower / upper): no data available
- h) Flash point: no data available
- i) Flashing temperature: no data available
- j) Decomposition temperature: no data available
- k) pH value: 0
- l) Kinematic viscosity: no data available
- m) Solubility in water: 0-100 %
- n) Dispersion coefficient ( $K_{o/w}$ ): no data available
- o) Vapour pressure (20°C): no data available
- p) Specific gravity: 1,11 g/cm<sup>3</sup>
- q) Relative vapour density ( $air=1$ ): no data available
- r) Particle size: no data available

**110 mL total Phosphate 15 Reagent A**

- a) State of aggregation: liquid
- b) Colour: colourless
- c) Odor: odorless
- d) Melting point: no data available
- e) Boiling point: no data available
- f) Flammability: no data available
- g) Explosive limits (lower / upper): no data available
- h) Flash point: no data available
- i) Flashing temperature: no data available
- j) Decomposition temperature: no data available
- k) pH value: 0-1
- l) Kinematic viscosity: no data available
- m) Solubility in water: 0-100 %
- n) Dispersion coefficient ( $K_{o/w}$ ): no data available
- o) Vapour pressure (20°C): no data available
- p) Specific gravity: 1,01 g/cm<sup>3</sup>
- q) Relative vapour density ( $air=1$ ): no data available
- r) Particle size: no data available

**2 g total Phosphate (B1)**

- a) State of aggregation: powder (solid)
- b) Colour: white
- c) Odor: odorless
- d) Melting point: no data available
- e) Boiling point: no data available
- f) Flammability: no data available
- g) Explosive limits (lower / upper): no data available
- h) Flash point: no data available
- i) Flashing temperature: no data available
- j) Decomposition temperature: no data available
- k) pH value: [10%] 6-7.5
- l) Kinematic viscosity: no data available
- m) Solubility in water: 0-38 %
- n) Dispersion coefficient ( $K_{o/w}$ ): no data available
- o) Vapour pressure (20°C): no data available
- p) Specific gravity: no data available
- q) Relative vapour density ( $air=1$ ): no data available



MACHEREY-NAGEL GmbH & Co. KG  
 Valenciennes Str. 11  
 52355 Düren · Germany  
[www.mn-net.com](http://www.mn-net.com)

DE Tel.: +49 24 21 969-0 info@mn-net.com  
 CH Tel.: +41 62 388 55 00 sales-ch@mn-net.com  
 FR Tel.: +33 388 68 22 68 sales-fr@mn-net.com  
 US Tel.: +1 888 321 62 24 sales-us@mn-net.com



# Safety Data Sheet

according to Regulations REACH 1907/2006/EC

REF: 985680

NANOCOLOR tot.Phosphate 15, robot

Page: 9/14

Printing date: 15.05.2024

Date of issue: 26.09.2022

Version: 2.2.3.10

r) Particle size: no data available

**5 mL total Phosphate (B2)**

a) State of aggregation:	liquid
b) Colour:	colourless
c) Odor:	odorless
d) Melting point:	no data available
e) Boiling point:	no data available
f) Flammability:	no data available
g) Explosive limits (lower / upper):	no data available
h) Flash point:	no data available
i) Flashing temperature:	no data available
j) Decomposition temperature:	no data available
k) pH value:	no data available
l) Kinematic viscosity:	no data available
m) Solubility in water:	no data available
n) Dispersion coefficient ( $K_{o/w}$ ):	no data available
o) Vapour pressure (20°C):	no data available
p) Specific gravity:	no data available
q) Relative vapour density ( $air=1$ ):	no data available
r) Particle size:	no data available

## 9.2 Other information

### 9.2.1 Information on physical hazard classes

no data available

### 9.2.2 Other safety-related parameters

No data is available for the other parameters for the mixtures, since no registration and no chemical safety report is required.

□ □

## SECTION 10: Stability and reactivity

### 10.1 Reactivity

no further data available.

### 10.2 Chemical stability

no known instability.

### 10.3 Possibility of hazardous reactions

No further data available.

### 10.4 Conditions to avoid

Persulfates decompose when heated by splitting off oxygen. Observe the storage temperature printed on it. No more required.

### 10.5 Incompatible materials

no additional data available

### 10.6 Hazardous decomposition products

In the original package all parts/all reagents are safety and separated stored. Decompositions are not observed during the expiration period under recommended conditions.

## SECTION 11: Toxicological information

### 11.1 Information on the hazard classes according regulation (EC) 1272/2008

Following information is valid for pure substances. Quantitative data on the toxicity of this product are not available.

**35 mg NANOFIX total Phosphate 15 (R0)**

Chemical:	<i>sodium peroxodisulfate</i>	CAS No.:	7775-27-1
TSCA Inventory:	listed	California Proposition 65 List:	not listed
Australia NICNAS:	Yes (PEC/18)	Canada CEPA 1999:	DSL Yes
Japan CSCL/PRTR:	not listed, Japan PDSCL:		not listed



MACHEREY-NAGEL GmbH & Co. KG  
 Valencienner Str. 11  
 52355 Düren · Germany  
[www.mn-net.com](http://www.mn-net.com)

DE Tel.: +49 24 21 969-0 info@mn-net.com  
 CH Tel.: +41 62 388 55 00 sales-ch@mn-net.com  
 FR Tel.: +33 388 68 22 68 sales-fr@mn-net.com  
 US Tel.: +1 888 321 62 24 sales-us@mn-net.com

# Safety Data Sheet

according to Regulations REACH 1907/2006/EC

REF: 985680	NANOCOLOR tot.Phosphate 15, robot	Page: 10/14
Printing date: 15.05.2024	Date of issue: 26.09.2022	Version: 2.2.3.10

Japan ISHL: listed  $\geq 1,0\%$ / $\geq 0,1\%$ , Article 57-2 (SDS required)  
 South Korea TCCA: not listed  
 Korea Exist.Chem.Inventory: KE-12369  
 LD50 orl rat : 902 mg/kg  
 Acute Effects: Cause after oral intake, inhalation of vapours/dust, skin contact, impairments of health when ingested in small quantities.  
 Chronic Effects: May cause sensitization by skin contact, also in repeated contact of small amounts. May cause allergy or asthma symptoms or breathing difficulties if inhaled.

**6 mL total Phosphate (C)**

Chemical: *ammonium heptamolybdate* CAS No.: 12054-85-2  
 TSCA Inventory: listed (CAS 11098-84-3)  
 Japan ISHL: listed  $\geq 1,0\%$ / $\geq 0,1\%$ ,  
 Korea Exist.Chem.Inventory: not listed  
 LD50 orl rat : 2000-5000 mg/kg  
 LD50 inh rat : 1,930-5,840 mg/L/4H

Chemical: *sulfuric acid* CAS No.: 7664-93-9  
 TSCA Inventory: listed California Proposition 65 List: not listed  
 ACGIH: 1 ppm  
 Exposure Routes: inhalation, ingestion, skin and/or eye contact  
 Target Organs: Eyes, skin, respiratory system, teeth  
 Symptoms: irritation eyes, skin, nose  
 Australia NICNAS: not listed Canada CEPA 1999: DSL Yes  
 Japan CSCL/PRTR: not listed, Japan PDSC: Deleterious Substance  
 Japan ISHL: listed  $\geq 1,0\%$ / $\geq 1,0\%$ , Article 57-2 (SDS required)  
 South Korea TCCA: Accident Precaution Chemical Yes  
 Korea Exist.Chem.Inventory: KE-32570  
 LD50 orl rat : 2140 mg/kg  
 LC50 inh mus : 0,85 mg/L/4H  
 TRGS 905 (DE): R F C

**110 mL total Phosphate 15 Reagent A**

Chemical: *sulfuric acid* CAS No.: 7664-93-9  
 TSCA Inventory: listed California Proposition 65 List: not listed  
 ACGIH: 1 ppm  
 Exposure Routes: inhalation, ingestion, skin and/or eye contact  
 Target Organs: Eyes, skin, respiratory system, teeth  
 Symptoms: irritation eyes, skin, nose  
 Australia NICNAS: not listed Canada CEPA 1999: DSL Yes  
 Japan CSCL/PRTR: not listed, Japan PDSC: Deleterious Substance  
 Japan ISHL: listed  $\geq 1,0\%$ / $\geq 1,0\%$ , Article 57-2 (SDS required)  
 South Korea TCCA: Accident Precaution Chemical Yes  
 Korea Exist.Chem.Inventory: KE-32570  
 LD50 orl rat : 2140 mg/kg  
 LC50 inh mus : 0,85 mg/L/4H  
 TRGS 905 (DE): R F C

**2 g total Phosphate (B1)**

Chemical: *sodium L(+)-ascorbate* CAS No.: 134-03-2  
 TSCA Inventory: listed  
 LD50 orl rat : 11900 mg/kg

**5 mL total Phosphate (B2)**

Chemical: *water* CAS No.: 7732-18-5  
 TSCA Inventory: listed  
 Korea Exist.Chem.Inventory: KE-35400  
 LD50 orl rat : > 90000 mg/kg



# Safety Data Sheet

according to Regulations REACH 1907/2006/EC

REF: 985680	NANOCOLOR tot.Phosphate 15, robot	Page: 11/14
Printing date: 15.05.2024	Date of issue: 26.09.2022	Version: 2.2.3.10

## 11.2 Other hazards

**Possible endocrine disrupting effects**  
no data available

**Other information**  
no additional data available

## SECTION 12: Ecological information

### 12.1 Toxicity

Following information is valid for pure substances.

#### 35 mg NANOFIX total Phosphate 15 (R0)

Substance name: *sodium peroxodisulfate* CAS-Nr.: 7775-27-1  
 Water hazard class (DE): 1 WGK No.: 1352  
 Storage class (VCI): 5.1 B

#### 6 mL total Phosphate (C)

Substance name: *ammonium heptamolybdate* CAS-Nr.: 12054-85-2  
 Water hazard class (DE): 1 WGK No.: 0637  
 Storage class (VCI): 12-13

Substance name: *sulfuric acid* CAS-Nr.: 7664-93-9

PNEC (fresh water): 2.5 µg/L  
 PNEC = Predicted No Effect Concentration = concentration at which no effect on the environment is expected

LC50 fish/96h : [NOEC, 65d] 25 µg/L  
 EC50 daphnia/48h : 100 mg/L  
 EC10 pseudomonas putida/16h : [72h] 100 mg/L  
 Water hazard class (DE): 1 WGK No.: 0182  
 Storage class (VCI): 8 B

#### 110 mL total Phosphate 15 Reagent A

Substance name: *sulfuric acid* CAS-Nr.: 7664-93-9

PNEC (fresh water): 2.5 µg/L  
 PNEC = Predicted No Effect Concentration = concentration at which no effect on the environment is expected

LC50 fish/96h : [NOEC, 65d] 25 µg/L  
 EC50 daphnia/48h : 100 mg/L  
 EC10 pseudomonas putida/16h : [72h] 100 mg/L  
 Water hazard class (DE): 1 WGK No.: 0182  
 Storage class (VCI): 8 B

#### 2 g total Phosphate (B1)

Substance name: *sodium L(+)-ascorbate* CAS-Nr.: 134-03-2

#### 5 mL total Phosphate (B2)

Substance name: *water* CAS-Nr.: 7732-18-5

### 12.2 Persistence and degradability

### 12.3 Bioaccumulative potential

### 12.4 Mobility in soil

### 12.5 Results of PBT and vPvB assessment

This substance/mixture contains no components considered to be either persistent, bioaccumulative and toxic (PBT) or very persistent and very bioaccumulative (vPvB) at levels of 0.1% or higher.

### 12.6 Endocrine disrupting properties

no data available



# Safety Data Sheet

according to Regulations REACH 1907/2006/EC

REF: 985680	NANOCOLOR tot.Phosphate 15, robot	Page: 12/14
Printing date: 15.05.2024	Date of issue: 26.09.2022	Version: 2.2.3.10

**12.7 Other adverse effects**  
no additional data available

## SECTION 13: Disposal considerations

Please observe local regulations for collection and disposal of hazardous waste and contact waste disposal company, where you will obtain information on laboratory waste disposal (waste code number 16 05 06).

**13.1 Waste treatment methods**  
Not necessary, see above.

## SECTION 14: Transport information

**14.1 - 14.4:** No dangerous goods according to the transport regulations because  
**De Minimis** (excepted quantities:  $\leq 1 \text{ mL} / \Sigma \leq 100 \text{ mL}$ ) = ADR 3.5.1.4  
**De Minimis** (excepted quantities:  $\leq 1 \text{ mL} / \Sigma \leq 100 \text{ mL}$ ) = IATA DRG 2.6.10

**14.5 Environmental hazards**  
none, contains only small quantities of hazardous substances

**14.6 Special precautions for user**  
not necessary

**14.7 Carriage in bulk by sea in accordance with IMO instruments**  
Not applicable.

## SECTION 15: Regulatory information

**15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**  
 Dangerous Substances Protection Act (DE: Chemikaliengesetz - ChemG), Aug 2013, Stand: Okt 2020  
 Ordinance on protection against dangerous substances (E: Gefahrstoffverordnung - GefStoffV), Nov 2010, Stand: Mrz 2017  
 TRGS 201, Classification and labeling of activities involving hazardous substances, Feb 2017  
 TRGS 220, National aspects when preparing safety data sheets, Jan 2017  
 TRGS 400, Risk assessment for activities involving hazardous substances, Jul 2017  
 TRGS 401, Skin contact hazard - identification, assessment, action, Jun 2008, status: Feb 2011  
 BekGS 408, Application of the GefStoffV and the TRGS with the entry into force of the CLP regulation, December 2009, status: Jan 2012  
 TRGS 500, Protective measures, Mai 2008  
 TRGS 510, Storage of hazardous substances in portable containers from March 2013, status: Oct 2015  
 Wasserhaushaltsgesetz - WHG, Section 3 Handling substances hazardous to water, Jul 2009, status: Aug 2016  
 MN leaflet/instructions for use, also at [www.mn-net.com](http://www.mn-net.com)  
 If necessary, observe other country-specific regulations.

**15.2 Chemical safety assessment**  
not necessary for these small amounts

## SECTION 16: Other information

**16.1 Changes compared to the last version**  
Between versions 2.2.3.10 and 2.2.2.2 following changes were applied: - 1 composition data corrected - 8 substance data corrected

### 16.2 List of H and P phrases

**16.2.1 List of relevant H phrases**

H	Between versions 2.2.3.10 and 2.2.2.2 following changes were applied: - 1 composition data corrected - 8 substance data corrected
H272	May intensify fire; oxidizer.
H302	Harmful if swallowed.
H315	Causes skin irritation.
H317	May cause an allergic skin reaction.
H319	Causes serious eye irritation.
H334	May cause allergy or asthma symptoms or breathing difficulties if inhaled.
H335	May cause respiratory irritation.

**16.2.2 List of relevant P phrases**



MACHEREY-NAGEL GmbH & Co. KG  
 Valencienner Str. 11  
 52355 Düren · Germany  
[www.mn-net.com](http://www.mn-net.com)

DE Tel.: +49 24 21 969-0 [info@mn-net.com](mailto:info@mn-net.com)  
 CH Tel.: +41 62 388 55 00 [sales-ch@mn-net.com](mailto:sales-ch@mn-net.com)  
 FR Tel.: +33 388 68 22 68 [sales-fr@mn-net.com](mailto:sales-fr@mn-net.com)  
 US Tel.: +1 888 321 62 24 [sales-us@mn-net.com](mailto:sales-us@mn-net.com)

# Safety Data Sheet

according to Regulations REACH 1907/2006/EC

REF: 985680	NANOCOLOR tot.Phosphate 15, robot	Page: 13/14
Printing date: 15.05.2024	Date of issue: 26.09.2022	Version: 2.2.3.10

P261sh	Avoid breathing dust/vapours.
P280sh	Wear protective gloves/eye protection.
P284	[In case of inadequate ventilation] wear respiratory protection.
P302+352	IF ON SKIN: Wash with plenty of water.
P333+313	If skin irritation or rash occurs: Get medical advice/attention.
P342+311	If experiencing respiratory symptoms: Call a POISON CENTER/doctor.
P362+364	Take off contaminated clothing and wash it before reuse.
P501	Dispose of contents/container to regulated waste treatment.

## 16.3 Recommended restriction on use

Only for professional user.  
 Look about employee restrictions for young people (f. ex. 94/33/EC or DE § 22 JArbSchG)!  
 Look about employee restrictions for pregnant women and nursing women (f.ex. 92/85/EEC or for DE §§ 11-13 MuSchG 2017)!  
 An individual package of this product or test kit has a moderate hazardous potential.

## 16.4 Sources of key data

KÜHN, BIRETT, Leaflets on hazardous materials, 2021  
 Directive 1999/92/EG Minimum requirements to improve the safety and health protection of workers at risk from potentially explosive atmospheres  
 SUVA .CH, limit values in the air at work 2009, revised on 01/2009  
 Regulation 790/2009/EU, adaptation of Regulation 1272/2008/EU to technical and scientific progress (1st ATP)  
 Regulation 453/2010/EU, adaptation of the REACH regulation 1907/2006/EG  
 TRGS 907, German technical rules for listing substances and causes of sensitization, updated November 2011  
 Regulation 487/2013/EU, adaptation of regulation 1272/2008/EG to technical and scientific progress (4th ATP)  
 Regulation 1221/2015/EU, adaptation of regulation 1272/2008/EG to technical and scientific progress (7th ATP)  
 Regulation 776/2017/EU, adaptation of regulation 1272/2008/EG to technical and scientific progress (10th ATP)  
  
 Regulation 669/2018/EU, adaptation of Regulation 1272/2008/EC to technical and scientific progressText (11th ATP)  
 Regulation 1480/2018/EU, adaptation of regulation 1272/2008/EG to technical and scientific progress (13th ATP)  
 Regulation 521/2019/EU, adaptation of regulation 1272/2008/EG to technical and scientific progress (12th ATP)  
 TRGS 900, German rules of technology on limit values in the air at work, as of 03/2019  
 Regulation 217/2020/EU, adaptation of Annex VI, Part 3, of Regulation 1272/2008/EC to technical and scientific progress (14th ATP)  
 Regulation 878/2020/EU, adaptation of Annex II of the REACH regulation 1907/2006/EG  
 Regulation 1182/2020/EU, adaptation of Annex VI, Part 3, of Regulation 1272/2008/EC to technical and scientific progress (15th ATP)  
 Regulation 643/2021/EU, adaptation of Annex VI, Part 1, of Regulation 1272/2008/EC to technical and scientific progress (16th ATP)  
 Regulation 849/2021/EU, adaptation of Annex VI, Part 3, of Regulation 1272/2008/EC to technical and scientific progress (17th ATP)  
 Regulation 692/2022/EU, adaptation of Annex VI, Part 1, of Regulation 1272/2008/EC to technical and scientific progress (18th ATP)

### revisions/updates

*Reason for revision:* 2014-02 Corrected structure of the sections according to Regulation 453/2010/EU, if necessary  
 2014-04 adjustment according Regulation 487/2013/EU  
 2016-03 adjustment according Regulation 1221/2015/EU  
  
 2017-11 adjustment according the ECHA registration dossier  
 2022-11 adjustment according Regulation 878/2020/EU

## 16.5 Further information

MACHEREY-NAGEL GmbH & Co. KG provides the information contained herein in good faith being up-to-date of own realizations at revision time. This document is intended only as a guide to the appropriate precautionary handling of the material by a properly trained person using this product. Individuals receiving the information must exercise their independent judgement in determining its appropriateness for a particular purpose.  
 MACHEREY-NAGEL GmbH & Co. KG makes NO REPRESENTATIONS or WARRANTIES, either expressed or implied, including without limitation any warranties of merchantability, fitness for a particular purpose with respect to the information set forth herein or the product to which the information refers. Accordingly MACHEREY-NAGEL GmbH & Co. KG will not be responsible for damages resulting from use of or reliance upon this information. See terms and conditions at the end of our price lists for additional information.

## 16.6 Legend / Abbreviations

acc:	according
ADR:	Convention concerning the International Carriage of Dangerous Goods by Road
Act:	acute
BAT:	biological workplace tolerance value
CAO:	Cargo Aircraft Only
Carc:	carcinogen
CAS:	Chemical Abstracts Service
CLP:	Classification, Labelling and Packaging regulation
CMR:	carcinogen, mutagen, reproduction toxic
Corr:	corrosive
COD:	chemical oxigen demand
CSSL:	Chemical Substance Control Law (Jp)



MACHEREY-NAGEL GmbH & Co. KG  
 Valenciener Str. 11  
 52355 Düren · Germany  
[www.mn-net.com](http://www.mn-net.com)

DE Tel.: +49 24 21 969-0 info@mn-net.com  
 CH Tel.: +41 62 388 55 00 sales-ch@mn-net.com  
 FR Tel.: +33 388 68 22 68 sales-fr@mn-net.com  
 US Tel.: +1 888 321 62 24 sales-us@mn-net.com

# Safety Data Sheet

according to Regulations REACH 1907/2006/EC

REF: 985680	NANOCOLOR tot.Phosphate 15, robot	Page: 14/14
Printing date: 15.05.2024	Date of issue: 26.09.2022	Version: 2.2.3.10

Dam:	damage
DNEL:	Derived No-Effect Level (for workers)
derm:	dermal
dog:	dog
EC10:	Concentration causing a toxic effect in 10% of the test organisms
EC:	European Community
EC-Nr:	Substance number of the EC substance inventory
EmS:	Guide to accident management measures on ships
EU:	European Union
fish:	fish (not specified)
GHS:	Global Harmonized System of Classification and Labeling of Chemicals
gpg:	guinea pig
ICAO:	International Civil Aviation Organization
ihl:	inhaled
IMDG:	International Maritime Dangerous Goods Code
intrav:	intravenous
ipt:	intraperitoneal
ISHL:	Industrial Safety and Health Law (Jp)
LC50:	letale concentration 50%
LD50:	letale dosis 50%
leuciscus idus:	fisch, ide, orfe
MAK:	maximum workplace concentration
Met:	Metall
mus:	mouse
Muta:	mutagen
NIOSH:	National Institute for Occupational Safety and Health (US)
NRD:	Non-rapidly degradable
onchorhynchus mykiss:	fish, rainbow trout
orl:	oral
OSHA:	Occupational Safety and Health Administration
PAX:	transport on passenger planes allowed
PBT:	persistent, bioaccumulating, toxic substance
pH:	pH value
pimephales promelas:	fish, fathead minnow
PNEC:	Predicted No Effect Concentration
PROC 15:	Process category 'for laboratory use'
PRTR:	Law for PRTR and Promotion of Chemical Management (Jp)
PVC:	polyvinyl chloride
quail:	bird, quail
rat:	rat
rbt:	rabbit
RD:	rapidly degradable
RE:	repeated
REACH:	Registration, Evaluation, Authorisation and Restriction of Chemicals
REF:	item number, reference number
Reg.No.:	rRegistration number
Repr:	harmful to reproduction
Resp:	respiratory
RIP:	REACH Implementations Projects
scu:	sub cutan
SDS:	safety data sheet
Sens:	sensitisation
STEL:	short term exposure limit
STOT:	Specific Target Organ Toxicity
SVHC:	Substance of Very High Concern
t/a:	tons per year
TCCA:	Toxic Chemicals Control Act (S. Korea)
Tox:	toxic
TSCA:	The Toxic Substances Control Act (US)
TWA:	time weighted average
TRGS:	technical regulations (DE)
vPvB:	very persistent, very bioaccumulating substance

## 16.7 Training advice

Regular safety training. Multiple safety training of staffs about danger and protection by using hazards in working area. Additionally training and introduction of staffs for using these products.



MACHEREY-NAGEL GmbH & Co. KG  
 Valenciener Str. 11  
 52355 Düren · Germany  
[www.mn-net.com](http://www.mn-net.com)

DE Tel.: +49 24 21 969-0 info@mn-net.com  
 CH Tel.: +41 62 388 55 00 sales-ch@mn-net.com  
 FR Tel.: +33 388 68 22 68 sales-fr@mn-net.com  
 US Tel.: +1 888 321 62 24 sales-us@mn-net.com