

DB-1

| ID (mm) | Length (m) | Film (µm) | Temp Limits (°C) | 7 in Cage | 5 in Cage | 7890/6890 | |
|---------|------------|----------------|------------------|-----------|-----------|---------------|-------------|
| | | | | | | LTM II Module | |
| 0.32 | 15 | 0.10 | -60 to 325/350 | 123-1011 | | 123-1011LTM | |
| | | 0.25 | -60 to 325/350 | 123-1012 | | 123-1012LTM | |
| | | 1.00 | -60 to 325/350 | 123-1013 | | | |
| | | 3.00 | -60 to 280/300 | 123-1014 | | | |
| | | 5.00 | -60 to 280/300 | 123-1015 | | 123-1015LTM | |
| | 25 | 0.12 | -60 to 325/350 | 123-1027 | | | |
| | | 0.25 | -60 to 325/350 | 123-1022 | | | |
| | | 0.52 | -60 to 325/350 | 123-1026 | | | |
| | | 1.05 | -60 to 325/350 | 123-102F | | | |
| | 30 | 0.10 | -60 to 325/350 | 123-1031 | | | |
| | | 0.25 | -60 to 325/350 | 123-1032 | | | 123-1032LTM |
| | | 0.50 | -60 to 325/350 | 123-103E | | | 123-103ELTM |
| | | 1.00 | -60 to 325/350 | 123-1033 | 123-1033E | | 123-1033LTM |
| | | 1.50 | -60 to 300/320 | 123-103B | | | 123-103BLTM |
| | | 3.00 | -60 to 280/300 | 123-1034 | | | |
| | | 5.00 | -60 to 280/300 | 123-1035 | | | 123-1035LTM |
| | 50 | 0.25 | -60 to 325/350 | 123-1052 | | | |
| | | 0.52 | -60 to 325/350 | 123-1056 | | | |
| | | 1.05 | -60 to 325/350 | 123-105F | | | |
| | | 1.20 | -60 to 325/350 | 123-105C | | | |
| 5.00 | | -60 to 280/300 | 123-1055 | | | | |
| 60 | 0.10 | -60 to 325/350 | 123-1061 | | | | |
| | 0.25 | -60 to 325/350 | 123-1062 | 123-1062E | | | |
| | 0.50 | -60 to 325/350 | 123-106E | | | | |
| | 1.00 | -60 to 325/350 | 123-1063 | 123-1063E | | | |
| | 1.50 | -60 to 300/320 | 123-106B | 123-106BE | | | |
| | 2.00 | -60 to 280/300 | 123-106G | | | | |
| | 3.00 | -60 to 280/300 | 123-1064 | 123-1064E | | | |
| | 5.00 | -60 to 280/300 | 123-1065 | 123-1065E | | | |
| 0.45 | 30 | 1.27 | -60 to 325/350 | 124-1032 | | | |
| | | 2.55 | -60 to 260/280 | 124-1034 | | | |

(Continued)