

**Agilent J&W LTM II Low Thermal Mass Column Modules
for 7890A/B Series GC Systems**

| Phase | ID (mm) | Length (m) | Film (μm) | Part No. |
|----------|---------|------------|------------------------|-------------|
| DB-1ht | 0.25 | 30 | 0.10 | 122-1131LTM |
| | 0.32 | 5 | 0.25 | 123-1102LTM |
| | | 15 | 0.10 | 123-1111LTM |
| DB-1ms | 0.10 | 10 | 0.10 | 127-0112LTM |
| | | 20 | 0.40 | 127-0123LTM |
| | 0.18 | 20 | 0.18 | 121-0122LTM |
| | 0.20 | 25 | 0.33 | 128-0122LTM |
| | 0.25 | 15 | 0.25 | 122-0112LTM |
| | | 30 | 0.25 | 122-0132LTM |
| DB-200 | 0.25 | 30 | 0.25 | 122-2032LTM |
| | | | 0.50 | 122-2033LTM |
| DB-210 | 0.53 | 30 | 1.00 | 125-0232LTM |
| DB-225 | 0.25 | 15 | 0.25 | 122-2212LTM |
| | | 30 | 0.25 | 122-2232LTM |
| DB-225ms | 0.25 | 15 | 0.25 | 122-2912LTM |
| | | 30 | 0.25 | 122-2932LTM |
| DB-23 | 0.25 | 30 | 0.25 | 122-2332LTM |
| DB-2887 | 0.53 | 10 | 3.00 | 125-2814LTM |
| DB-35 | 0.32 | 30 | 0.50 | 123-1933LTM |
| | | | 1.00 | 125-1932LTM |
| DB-35ms | 0.25 | 30 | 0.25 | 122-3832LTM |
| DB-5 | 0.10 | 10 | 0.10 | 127-5012LTM |
| | | | 0.17 | 127-501ELTM |
| | | | 0.40 | 127-5013LTM |
| | 0.15 | 10 | 1.20 | 12A-5015LTM |
| | | | 0.18 | 121-5012LTM |
| | 0.18 | 10 | 0.40 | 121-5013LTM |
| | | | 0.18 | 121-5022LTM |
| | | | 0.40 | 121-5023LTM |
| | 0.20 | 25 | 0.33 | 128-5022LTM |

(Continued)