

DB-5

ID (mm)	Length (m)	Film (μm)	Temp Limits ($^{\circ}\text{C}$)	7 in Cage	5 in Cage	7890/6890		
						LTM II Module		
0.45	30	0.42	-60 to 300/320	124-5037				
		1.27	-60 to 300/320	124-5032				
0.53	10	2.65	-60 to 260/280	125-50HB				
		15	0.25	-60 to 300/320	125-501K			
			0.50	-60 to 300/320	125-5017			
	1.00		-60 to 300/320	125-501J				
	25	5.00	1.50	-60 to 300/320	125-5012	125-5012E	125-5012LTM	
				-60 to 260/280	125-5025			
	30	0.25		-60 to 300/320	125-503K			
				-60 to 300/320	125-5037			
				-60 to 300/320	125-503D			
				-60 to 300/320	125-503J			
		1.50		-60 to 300/320	125-5032	125-5032E	125-5032LTM	
				-60 to 260/280	125-503B			
				-60 to 260/280	125-5034			
			-60 to 260/280	125-5035	125-5035E	125-5035LTM		
60	1.50	-60 to 300/320	125-5062					
	5.00	-60 to 260/280	125-5065	125-5065E				