

DB-WAXetr

- Polyethylene glycol (PEG)
- Extended temperature range (etr)
- High polarity
- Excellent column-to-column repeatability
- Bonded and cross-linked
- Solvent rinsable
- Equivalent to USP Phase G16

Similar Phases: SUPELCOWAX 10, SUPEROX II, CB-WAX, Stabilwax, BP-20, 007-CW, Carbowax, Rtx-WAX, ZB-WAX, ZB-WAX plus

DB-WAXetr

| ID (mm) | Length (m) | Film (μm) | Temp Limits ($^{\circ}\text{C}$) | 7 in Cage | 5 in Cage | 7890/6890 |
|---------|------------|------------------------|------------------------------------|-----------|-----------|---------------|
| | | | | | | LTM II Module |
| 0.20 | 25 | 0.40 | 30 to 250/260 | 128-7323 | | |
| 0.25 | 30 | 0.25 | 30 to 260/280 | 122-7332 | 122-7332E | 122-7332LTM |
| | | 0.50 | 30 to 250/260 | 122-7333 | | |
| | 60 | 0.25 | 30 to 260/280 | 122-7362 | | |
| | | 0.50 | 30 to 250/260 | 122-7363 | | |
| 0.32 | 15 | 0.25 | 30 to 260/280 | 123-7312 | | |
| | | 1.00 | 30 to 250/260 | 123-7314 | | |
| | 30 | 0.25 | 30 to 260/280 | 123-7332 | | |
| | | 0.50 | 30 to 250/260 | 123-7333 | | |
| | | 1.00 | 30 to 250/260 | 123-7334 | | 123-7334LTM |
| | 50 | 1.00 | 30 to 250/260 | 123-7354 | 123-7354E | |
| | 60 | 0.25 | 30 to 260/280 | 123-7362 | | |
| | | 0.50 | 30 to 250/260 | 123-7363 | | |
| 1.00 | | 30 to 250/260 | 123-7364 | | | |
| 0.53 | 15 | 1.00 | 30 to 240/260 | 125-7312 | | |
| | | 2.00 | 50 to 230/250 | 125-7314 | | |
| | 30 | 1.00 | 30 to 240/260 | 125-7332 | 125-7332E | |
| | | 1.50 | 30 to 230/240 | 125-7333 | | 125-7333LTM |
| | | 2.00 | 50 to 230/250 | 125-7334 | 125-7334E | |
| | 60 | 1.00 | 30 to 240/260 | 125-7362 | | |