



Edelstahl-Fittings (S),  
5062-2418



PEEK-Fittings (SPF),  
0100-1516/5063-6591



Edelstahl-Fittings, lang (SL),  
5065-4454



PEEK-Fitting (SPF), fingerfest,  
0100-1516



Edelstahl-Fitting, extralang (SX),  
5065-9967



Lange PEEK-Fittings (SPFL),  
5062-8541



Edelstahl-Fitting,  
5061-3303



PEEK-Fittings (SPF), fingerfest,  
5065-4426



Vordere Ferrulen, Edelstahl,  
5180-4108



Doppelflügel-Fitting (SPF),  
5042-6500



Hintere Ferrule, Edelstahl,  
5180-4114



PEEK RheFlex-Fittings (SPF),  
0100-1631



Austauschbares Fitting (SV), 1200 bar,  
5067-4733



PEEK RheFlex-Fittings (SPF),  
0100-2175



Langes austauschbares Fitting (SLV),  
1200 bar, 5067-4738



Verschlusskappe, Edelstahl,  
01080-83202



Extralanges austauschbares Fitting  
(SXV), 1200 bar, 5067-4739

## Fittings

| Beschreibung                                                            | Legende | Einheit | Best.-Nr.   |
|-------------------------------------------------------------------------|---------|---------|-------------|
| Swagelok-Edelstahl-Fitting, 1,6 mm                                      | S       | 10 St.  | 5062-2418   |
| Swagelok-Edelstahl-Fitting, 1,6 mm, lange Schraube                      | SL      | 10 St.  | 5065-4454   |
| Swagelok-Edelstahl-Fitting, 1,6 mm, extralange Schraube                 | SX      | 10 St.  | 5065-9967   |
| Swagelok-1,6 mm-Schraube                                                |         | 10 St.  | 5061-3303   |
| Vordere Ferrule, 1,6 mm, Edelstahl                                      |         | 10 St.  | 5180-4108   |
| Hintere Ferrule, 1,6 mm, Edelstahl                                      |         | 10 St.  | 5180-4114   |
| Austauschbares Swagelok-Fitting, 1,6 mm, 1200 bar                       | SV      |         | 5067-4733   |
| Austauschbares Swagelok-Fitting, 1,6 mm, 1200 bar, lange Schraube       | SLV     |         | 5067-4738   |
| Austauschbares Swagelok-Fitting, 1,6 mm, 1200 bar, extralange Schraube  | SXV     |         | 5067-4739   |
| PEEK Swagelok-Fitting, 1,6 mm, fingerfest                               | SPF     | 10 St.  | 5063-6591   |
| PEEK Swagelok-Fitting, 1,6 mm, fingerfest                               | SPF     | 2 St.   | 0100-1516   |
| Langes PEEK Swagelok-Fitting, 1,6 mm, fingerfest                        | SPLF    | 10 St.  | 5062-8541   |
| PEEK Swagelok-Fitting, 1,6 mm, fingerfest (verschiedene Farben)         | SPF     | 10 St.  | 5065-4426   |
| PEEK Swagelok-Doppelflügel-Fitting, 1,6 mm, fingerfest                  | SPF     | 10 St.  | 5042-6500   |
| PEEK Swagelok RheFlex-Fitting, 1,6 mm, fingerfest                       | SPF     | 5 St.   | 0100-1631   |
| PEEK Swagelok RheFlex-Fitting, 1,6 mm, fingerfest (verschiedene Farben) | SPF     | 10 St.  | 0100-2175   |
| Swagelok-Verschlusskappe, 1,6 mm, Edelstahl                             | S       | je 1    | 01080-83202 |

## Material

| Legende | Beschreibung                                          |
|---------|-------------------------------------------------------|
| S       | Swagelok, 1,6 mm Anschluss-ID                         |
| SL      | Swagelok, 1,6 mm Anschluss-ID, lang                   |
| SX      | Swagelok, 1,6 mm Anschluss-ID, extralang              |
| SV      | Swagelok, 1,6 mm Anschluss-ID, 1200 bar               |
| SLV     | Swagelok, 1,6 mm Anschluss-ID, lang, 1200 bar         |
| SLB     | Swagelok, 1,6 mm Anschluss-ID, lang, bio              |
| SXV     | Swagelok, 1,6 mm Anschluss-ID, extralang, 1200 bar    |
| SPF     | Swagelok, 1,6 mm Anschluss-ID, PEEK, fingerfest       |
| SPLF    | Swagelok, 1,6 mm Anschluss-ID, PEEK, lang, fingerfest |