



Edelstahl-Fittings (S),  
5062-2418



PEEK-Fittings (SPF),  
0100-1516/5063-6591



Edelstahl-Fittings, lang (SL),  
5065-4454



PEEK-Fitting (SPF), fingerfest,  
0100-1516



Edelstahl-Fitting, extralange (SX),  
5065-9967



Lange PEEK-Fittings (SPFL),  
5062-8541



Edelstahl-Fitting,  
5061-3303



PEEK-Fittings (SPF), fingerfest,  
5065-4426



Vordere Ferrulen, Edelstahl,  
5180-4108



Doppelflügel-Fitting (SPF),  
5042-6500



Hintere Ferrule, Edelstahl,  
5180-4114



PEEK RheFlex-Fittings (SPF),  
0100-1631



Austauschbares Fitting (SV), 1200 bar,  
5067-4733



PEEK RheFlex-Fittings (SPF),  
0100-2175



Langes austauschbares Fitting (SLV),  
1200 bar, 5067-4738



Verschlusskappe, Edelstahl,  
01080-83202



Extralanges austauschbares Fitting  
(SXV), 1200 bar, 5067-4739

## Fittings

Beschreibung	Legende	Einheit	Best.-Nr.
Swagelok-Edelstahl-Fitting, 1,6 mm	S	10 St.	5062-2418
Swagelok-Edelstahl-Fitting, 1,6 mm, lange Schraube	SL	10 St.	5065-4454
Swagelok-Edelstahl-Fitting, 1,6 mm, extralange Schraube	SX	10 St.	5065-9967
Swagelok-1,6 mm-Schraube		10 St.	5061-3303
Vordere Ferrule, 1,6 mm, Edelstahl		10 St.	5180-4108
Hintere Ferrule, 1,6 mm, Edelstahl		10 St.	5180-4114
Austauschbares Swagelok-Fitting, 1,6 mm, 1200 bar	SV		5067-4733
Austauschbares Swagelok-Fitting, 1,6 mm, 1200 bar, lange Schraube	SLV		5067-4738
Austauschbares Swagelok-Fitting, 1,6 mm, 1200 bar, extralange Schraube	SXV		5067-4739
PEEK Swagelok-Fitting, 1,6 mm, fingerfest	SPF	10 St.	5063-6591
PEEK Swagelok-Fitting, 1,6 mm, fingerfest	SPF	2 St.	0100-1516
Langes PEEK Swagelok-Fitting, 1,6 mm, fingerfest	SPLF	10 St.	5062-8541
PEEK Swagelok-Fitting, 1,6 mm, fingerfest (verschiedene Farben)	SPF	10 St.	5065-4426
PEEK Swagelok-Doppelflügel-Fitting, 1,6 mm, fingerfest	SPF	10 St.	5042-6500
PEEK Swagelok RheFlex-Fitting, 1,6 mm, fingerfest	SPF	5 St.	0100-1631
PEEK Swagelok RheFlex-Fitting, 1,6 mm, fingerfest (verschiedene Farben)	SPF	10 St.	0100-2175
Swagelok-Verschlusskappe, 1,6 mm, Edelstahl	S	je 1	01080-83202

## Material

Legende	Beschreibung
S	Swagelok, 1,6 mm Anschluss-ID
SL	Swagelok, 1,6 mm Anschluss-ID, lang
SX	Swagelok, 1,6 mm Anschluss-ID, extralange
SV	Swagelok, 1,6 mm Anschluss-ID, 1200 bar
SLV	Swagelok, 1,6 mm Anschluss-ID, lang, 1200 bar
SLB	Swagelok, 1,6 mm Anschluss-ID, lang, bio
SXV	Swagelok, 1,6 mm Anschluss-ID, extralange, 1200 bar
SPF	Swagelok, 1,6 mm Anschluss-ID, PEEK, fingerfest
SPLF	Swagelok, 1,6 mm Anschluss-ID, PEEK, lang, fingerfest