

Hi-Plex-Säulen zur Kohlenhydratanalyse

Beschreibung	LC Säulenabm. (mm)	Partikelgröße (µm)	Quervernetzungen (%)	Gegenion	Best.-Nr.
Hi-Plex Ca USP L19	4,0 x 250	8	8	Ca ²⁺	PL1570-5810
Hi-Plex Ca (Duo)	6,5 x 300	8	8	Ca ²⁺	PL1F70-6850
Hi-Plex Ca	7,7 x 300	8	8	Ca ²⁺	PL1170-6810
Hi-Plex Pb USP L34	7,7 x 100	8	8	Pb ²⁺	PL1170-2820
Hi-Plex Pb	7,7 x 300	8	8	Pb ²⁺	PL1170-6820
Hi-Plex K	7,7 x 300	8	8	K ⁺	PL1170-6860
Hi-Plex H	6,5 x 300	8	8	H ⁺	PL1F70-6830
Hi-Plex H	7,7 x 300	8	8	H ⁺	PL1170-6830
Hi-Plex H USP L17	7,7 x 100	8	8	H ⁺	PL1170-2823
Hi-Plex Na	7,7 x 300	10	4	Na ⁺	PL1171-6140
Hi-Plex Na (Octo)	7,7 x 300	8	8	Na ⁺	PL1170-6840

Hi-Plex Vorsäulen

Beschreibung	LC Säulenabm. (mm)	Partikelgröße (µm)	Quervernetzungen (%)	Gegenion	Best.-Nr.
Hi-Plex Ca	7,7 x 50	8	8	Ca ²⁺	PL1170-1810
Hi-Plex Ca (Duo)	7,7 x 50	8	8	Ca ²⁺	PL1170-1850
Hi-Plex Pb	7,7 x 50	8	8	Pb ²⁺	PL1170-1820
Hi-Plex K	7,7 x 50	8	8	K ⁺	PL1170-1860
Hi-Plex H	7,7 x 50	8	8	H ⁺	PL1170-1830
Hi-Plex Na	7,7 x 50	10	4	Na ⁺	PL1171-1140
Hi-Plex Na (Octo)	7,5 x 50	8	8	Na ⁺	PL1170-1840

Hi-Plex Vorsäulenkartuschen, 2 Stück

Beschreibung	LC Säulenabm. (mm)	Partikelgröße (µm)	Quervernetzungen (%)	Gegenion	Best.-Nr.
Hi-Plex Ca	3,0 x 5,0	8	8	Ca ²⁺	PL1670-0810
Hi-Plex Ca (Duo)	3,0 x 5,0	8	8	Ca ²⁺	PL1670-0850
Hi-Plex Pb	3,0 x 5,0	8	8	Pb ²⁺	PL1670-0820
Hi-Plex K	3,0 x 5,0	8	8	K ⁺	PL1670-0860
Hi-Plex H	3,0 x 5,0	8	8	H ⁺	PL1670-0830
Hi-Plex Na	3,0 x 5,0	10	4	Na ⁺	PL1671-0140
Hi-Plex Na (Octo)	3,0 x 5,0	8	8	Na ⁺	PL1670-0840
Guard Kartuschenhalter für 3,0 x 5,0 mm-Kartuschen					PL1310-0016