

Zubehör Accessories

7835



Sicherheits-Gasschlauch DIN 30664/ Schlauchselle (7760)

Längenempfehlung: 1 m, Preis pro Meter
Innendurchmesser x Wandstärke: 10 x 2 mm

Safety gas tube DIN 30664/ Tubing clamp (7760)

Recommendation of length: 1 m, price per meter
Inner diameter x wall thickness: 10 x 2 mm



Material	Artikel Nr.
Kautschuk/cautchouc	7835
Stahl vernickelt	7760


Spiritusbrenner 18/10 Stahl

Inhalt 150 ml, regulierbar, inkl. Docht

Alcohol burner 18/10 stainless steel

Volume 150 ml, adjustable, incl. wick

7200



Material	D	C	Artikel Nr.
18/10-Stahl	80	65	7200

Bedienungsanleitung als PDF zum Download auf www.bochem.de
Download of operating manual as PDF on www.bochem.com



7201

Docht für Spiritusbrenner

100 % Baumwolle

Wick for alcohol burner

100 % cotton

Material	A	Artikel Nr.
Wolle/Cotton	1 m	7201

Drahtdreiecke, Stahl verzinkt

Mit Tonrohren

Triangles, steel zincked

With clay tubes

Material	A	Artikel Nr.
Stahl verzinkt	50	12701
Stahl verzinkt	60	12702
Stahl verzinkt	70	12703

12701

