

## Tiegel und Deckel, Nickel 99,5 %

Deckel bitte separat bestellen.

### Crucible and lid, nickel 99,5 %

Please order lids separately.

Material	C	D	ML	Typ	Artikel Nr.
<b>Tiegel/Crucibles</b>					
Nickel 99,5 %	20	20	5	1	8100
Nickel 99,5 %	25	25	10	1	8101
Nickel 99,5 %	30	30	15	1	8102
Nickel 99,5 %	35	35	25	1	8103
Nickel 99,5 %	40	40	30	1	8104
Nickel 99,5 %	45	45	50	1	8105
Nickel 99,5 %	50	50	70	1	8106
Nickel 99,5 %	60	60	130	1	8107
Nickel 99,5 %	80	80	270	1	8109
Nickel 99,5 %	30	30	15	2	8122
Nickel 99,5 %	35	35	25	2	8123
Nickel 99,5 %	40	40	30	2	8124
Nickel 99,5 %	45	45	50	2	8125
Nickel 99,5 %	50	50	70	2	8126
Nickel 99,5 %	60	60	130	2	8127
Nickel 99,5 %	80	80	270	2	8129
Nickel 99,5 %	45	45	50	3	8145
Nickel 99,5 %	40	40	30	4	8164
Nickel 99,5 %	50	50	70	4	8166
Nickel 99,5 %	60	60	130	4	8167
<b>Deckel/Lids</b>					
Nickel 99,5 %		20			8180
Nickel 99,5 %		25			8181
Nickel 99,5 %		30			8182
Nickel 99,5 %		35			8183
Nickel 99,5 %		40			8184
Nickel 99,5 %		45			8185
Nickel 99,5 %		50			8186
Nickel 99,5 %		60			8187
Nickel 99,5 %		80			8189

## Tiegel und Deckel, 18/10-Stahl

Deckel bitte separat bestellen.

### Crucible and lid, 18/10 stainless steel

Please order lids separately.


Material	C	D	ML	Typ	Artikel Nr.
<b>Tiegel/Crucibles</b>					
18/10-Stahl	35	35	25	2	8881
18/10-Stahl	35	40	30	2	8882
18/10-Stahl	45	50	70	2	8883
<b>Deckel/Lids</b>					
18/10-Stahl		35			8891
18/10-Stahl		40			8892
18/10-Stahl		50			8893

Typ 1: 0,5 mm

Typ 2: 1,0 mm

Typ 3: 1,5 mm

Typ 4: 2,0 mm



2 verschieden Materialien:

Nickel 99,5 %

18/10-Stahl

2 different materials:

Nickel 99,5 %

18/10 stainless steel

