

Typ 1



5423

Typ 2



5553


Stativklemme 4 Finger, 18/10-Stahl

Finger mit Silikonüberzug


Retort clamp 4 prongs, 18/10 stainless steel

Finger with silicone-coating

Material	A	D	d	Typ	Artikel Nr.
18/10-Stahl	150	12	0-120	1	5423
18/10-Stahl	190	12	30-100	2	5553
18/10-Stahl	190	12	50-150	2	5554



B08297100



Stativklemme 4 Finger, 18/10-Stahl/Silikon

Spannweite d durch 2 Stellschrauben unabhängig voneinander regelbar.

Retort clamp 4 prong, 18/10 steel/silicone

Span width d independently adjustable by 2 regulating screws.

Material	A	D	d	Artikel Nr.
18/10-Stahl/Silikon	140	12	35-100	B08297100

Typ 1



5400

Typ 2



5405

Mikroklemme, 18/10-Stahl

Typ 1+2: Finger mit Kork

Typ 3: Finger mit Silikon

Clamp micro, 18/10 stainless steel

Type 1+2: Finger with cork

Type 3: Finger with silicone

Material	A	D	d	Typ	Artikel Nr.
18/10-Stahl	150	10	15	1	5400
18/10-Stahl	150	10	15	2	5405
18/10-Stahl	150	10	0-20	3	5415

Typ 3



5415

