

## Stativklemme Standard

\*Spannweite d durch 2 Stellschrauben unabhängig voneinander regelbar.

- Laboral/PVC
- 18/10-Stahl/Silikon

## Retort clamp standard

\*Span width d independently adjustable by 2 regulating screws.


- Laboral/PVC
- 18/10 steel/silicone

Material	A	D	d	Artikel Nr.
Laboral/PVC	125	8	8-30	B08235037
18/10-Stahl/Silikon	125	8	8-30	B08307036*
Laboral/PVC	140	12	15-65	B08257658*
18/10-Stahl/Silikon	140	12	20-75	B08294075
Laboral/PVC	140	12	22-70	B08301071

PVC = Polyvinylchlorid/polyvinyl chloride, -20/+80 °C



B08235036



5703



5704

## 2-Finger-Mehrzweckklemme, 18/10-Stahl NEW!

Spannweite d durch 2 Schrauben unabhängig voneinander regelbar, bietet eine große Spannweite.

Finger mit rutschfestem Vinylüberzug versehen.

## 2-Prong Multi-Purpose Clamps, 18/10 stainless steel

Span width d independently adjustable by 2 regulating screws, offers a wide range of motion. Prongs supplied with non-slip vinyl sleeves.

Material	A	D	d	Artikel Nr.
18/10-Stahl	127	11	0-75	5703
18/10-Stahl	127	11	0-95	5704



5551



5420

## Stativklemme 3 Finger, 18/10-Stahl

Finger mit Silikonüberzug

## Retort clamp 3 prongs, 18/10 stainless steel

Finger with silicone-coating

Material	A	D	d	Artikel Nr.
18/10-Stahl	190	12	0-80	5551
18/10-Stahl	190	12	0-120	5420

