

| REFERENCE | PHASE | LENGTH (mm) | ID (mm) | PARTICLE SIZE (µ) | TYPE |
|-----------|-------------------|-------------|---------|-------------------|-----------------|
| 20042 | CHIRALPAK® AS | 50 | 20 | 10 | Guard Semi-Prep |
| 20045 | CHIRALPAK® AS | 250 | 20 | 10 | Semi-Prep |
| 20094 | CHIRALPAK® AS | 150 | 2.1 | 10 | Analytical |
| 20135 | CHIRALPAK® AS SFC | 250 | 10 | 10 | Semi-Prep |
| 20145 | CHIRALPAK® AS SFC | 250 | 20 | 10 | Semi-Prep |
| 20175 | CHIRALPAK® AS SFC | 250 | 30 | 10 | Prep |

CHIRALPAK® AS-H

| | | | | | |
|-------|---------------------|-----|-----|---|------------------------|
| 20311 | CHIRALPAK® AS-H | 10 | 4 | 5 | Guard Cartridges (x 3) |
| 20321 | CHIRALPAK® AS-H | 50 | 4.6 | 5 | Analytical |
| 20323 | CHIRALPAK® AS-H | 100 | 4.6 | 5 | Analytical |
| 20324 | CHIRALPAK® AS-H | 150 | 4.6 | 5 | Analytical |
| 20325 | CHIRALPAK® AS-H | 250 | 4.6 | 5 | Analytical |
| 20335 | CHIRALPAK® AS-H | 250 | 10 | 5 | Semi-Prep |
| 20337 | CHIRALPAK® AS-H | 20 | 10 | 5 | Guard Semi-Prep |
| 20345 | CHIRALPAK® AS-H | 250 | 20 | 5 | Semi-Prep |
| 20375 | CHIRALPAK® AS-H | 250 | 30 | 5 | Prep |
| 20355 | CHIRALPAK® AS-H | 250 | 50 | 5 | Prep |
| 20394 | CHIRALPAK® AS-H | 150 | 2.1 | 5 | Analytical |
| 20423 | CHIRALPAK® AS-H SFC | 100 | 4.6 | 5 | Analytical |
| 20435 | CHIRALPAK® AS-H SFC | 250 | 10 | 5 | Semi-Prep |
| 20445 | CHIRALPAK® AS-H SFC | 250 | 20 | 5 | Semi-Prep |
| 20475 | CHIRALPAK® AS-H SFC | 250 | 30 | 5 | Prep |
| 20455 | CHIRALPAK® AS-H SFC | 250 | 50 | 5 | Prep |

CHIRALPAK® AS-3

| | | | | | |
|-------|---------------------|-----|-----|---|------------------------|
| 20511 | CHIRALPAK® AS-3 | 10 | 4 | 3 | Guard Cartridges (x 3) |
| 20522 | CHIRALPAK® AS-3 | 50 | 4.6 | 3 | Analytical |
| 20523 | CHIRALPAK® AS-3 | 100 | 4.6 | 3 | Analytical |
| 20524 | CHIRALPAK® AS-3 | 150 | 4.6 | 3 | Analytical |
| 20525 | CHIRALPAK® AS-3 | 250 | 4.6 | 3 | Analytical |
| 20592 | CHIRALPAK® AS-3 | 50 | 2.1 | 3 | Analytical |
| 20593 | CHIRALPAK® AS-3 | 100 | 2.1 | 3 | Analytical |
| 20594 | CHIRALPAK® AS-3 | 150 | 2.1 | 3 | Analytical |
| 20595 | CHIRALPAK® AS-3 | 250 | 2.1 | 3 | Analytical |
| 20S82 | CHIRALPAK® AS-3 SFC | 50 | 3 | 3 | Analytical |
| 20S83 | CHIRALPAK® AS-3 SFC | 100 | 3 | 3 | Analytical |
| 20S84 | CHIRALPAK® AS-3 SFC | 150 | 3 | 3 | Analytical |

CHIRALPAK® AS-RH

| | | | | | |
|-------|------------------|-----|-----|---|------------------------|
| 20711 | CHIRALPAK® AS-RH | 10 | 4 | 5 | Guard Cartridges (x 3) |
| 20724 | CHIRALPAK® AS-RH | 150 | 4.6 | 5 | Analytical |
| 20794 | CHIRALPAK® AS-RH | 150 | 2.1 | 5 | Analytical |

CHIRALPAK® AS-3R

| | | | | | |
|-------|------------------|-----|-----|---|------------------------|
| 20811 | CHIRALPAK® AS-3R | 10 | 4 | 3 | Guard Cartridges (x 3) |
| 20822 | CHIRALPAK® AS-3R | 50 | 4.6 | 3 | Analytical |
| 20823 | CHIRALPAK® AS-3R | 100 | 4.6 | 3 | Analytical |
| 20824 | CHIRALPAK® AS-3R | 150 | 4.6 | 3 | Analytical |
| 20892 | CHIRALPAK® AS-3R | 50 | 2.1 | 3 | Analytical |
| 20893 | CHIRALPAK® AS-3R | 100 | 2.1 | 3 | Analytical |
| 20894 | CHIRALPAK® AS-3R | 150 | 2.1 | 3 | Analytical |
| 20895 | CHIRALPAK® AS-3R | 250 | 2.1 | 3 | Analytical |

| REFERENCE | PHASE | LENGTH (mm) | ID (mm) | PARTICLE SIZE (µ) | TYPE |
|-----------------------------|-------------------|-------------|---------|-------------------|------------|
| CHIRALPAK® AY (T304) | | | | | |
| 47020 | CHIRALPAK® AY | - | Kg | 20 | CSP |
| 47021 | CHIRALPAK® AY | - | 100 g | 20 | CSP |
| 47224 | CHIRALPAK® AY | 150 | 4.6 | 20 | Analytical |
| 47225 | CHIRALPAK® AY | 250 | 4.6 | 20 | Analytical |
| 47226 | CHIRALPAK® AY | 300 | 4.6 | 20 | Analytical |
| 47010 | CHIRALPAK® AY | - | Kg | 10 | CSP |
| 47024 | CHIRALPAK® AY | 150 | 4.6 | 10 | Analytical |
| 47025 | CHIRALPAK® AY | 250 | 4.6 | 10 | Analytical |
| 47135 | CHIRALPAK® AY SFC | 250 | 10 | 10 | Semi-Prep |
| 47145 | CHIRALPAK® AY SFC | 250 | 20 | 10 | Semi-Prep |
| 47175 | CHIRALPAK® AY SFC | 250 | 30 | 10 | Prep |

CHIRALPAK® AY-H

| | | | | | |
|-------|---------------------|-----|-----|---|------------------------|
| 47311 | CHIRALPAK® AY-H | 10 | 4 | 5 | Guard Cartridges (x 3) |
| 47321 | CHIRALPAK® AY-H | 50 | 4.6 | 5 | Analytical |
| 47323 | CHIRALPAK® AY-H | 100 | 4.6 | 5 | Analytical |
| 47324 | CHIRALPAK® AY-H | 150 | 4.6 | 5 | Analytical |
| 47325 | CHIRALPAK® AY-H | 250 | 4.6 | 5 | Analytical |
| 47335 | CHIRALPAK® AY-H | 250 | 10 | 5 | Semi-Prep |
| 47337 | CHIRALPAK® AY-H | 20 | 10 | 5 | Guard Semi-Prep |
| 47345 | CHIRALPAK® AY-H | 250 | 20 | 5 | Semi-Prep |
| 47375 | CHIRALPAK® AY-H | 250 | 30 | 5 | Prep |
| 47355 | CHIRALPAK® AY-H | 250 | 50 | 5 | Prep |
| 47394 | CHIRALPAK® AY-H | 150 | 2.1 | 5 | Analytical |
| 47423 | CHIRALPAK® AY-H SFC | 100 | 4.6 | 5 | Analytical |
| 47435 | CHIRALPAK® AY-H SFC | 250 | 10 | 5 | Semi-Prep |
| 47445 | CHIRALPAK® AY-H SFC | 250 | 20 | 5 | Semi-Prep |
| 47455 | CHIRALPAK® AY-H SFC | 250 | 50 | 5 | Prep |
| 47475 | CHIRALPAK® AY-H SFC | 250 | 30 | 5 | Prep |

CHIRALPAK® AY-3

| | | | | | |
|-------|---------------------|-----|-----|---|------------------------|
| 47511 | CHIRALPAK® AY-3 | 10 | 4 | 3 | Guard Cartridges (x 3) |
| 47522 | CHIRALPAK® AY-3 | 50 | 4.6 | 3 | Analytical |
| 47523 | CHIRALPAK® AY-3 | 100 | 4.6 | 3 | Analytical |
| 47524 | CHIRALPAK® AY-3 | 150 | 4.6 | 3 | Analytical |
| 47525 | CHIRALPAK® AY-3 | 250 | 4.6 | 3 | Analytical |
| 47592 | CHIRALPAK® AY-3 | 50 | 2.1 | 3 | Analytical |
| 47593 | CHIRALPAK® AY-3 | 100 | 2.1 | 3 | Analytical |
| 47594 | CHIRALPAK® AY-3 | 150 | 2.1 | 3 | Analytical |
| 47594 | CHIRALPAK® AY-3 | 250 | 2.1 | 3 | Analytical |
| 47S82 | CHIRALPAK® AY-3 SFC | 50 | 3 | 3 | Analytical |
| 47S83 | CHIRALPAK® AY-3 SFC | 100 | 3 | 3 | Analytical |
| 47S84 | CHIRALPAK® AY-3 SFC | 150 | 3 | 3 | Analytical |

CHIRALPAK® AY-RH

| | | | | | |
|-------|------------------|-----|-----|---|------------------------|
| 47711 | CHIRALPAK® AY-RH | 10 | 4 | 5 | Guard Cartridges (x 3) |
| 47724 | CHIRALPAK® AY-RH | 150 | 4.6 | 5 | Analytical |
| 47794 | CHIRALPAK® AY-RH | 150 | 2.1 | 5 | Analytical |

CHIRALPAK® AY-3R

| | | | | | |
|-------|------------------|----|-----|---|------------------------|
| 47811 | CHIRALPAK® AY-3R | 10 | 4 | 3 | Guard Cartridges (x 3) |
| 47822 | CHIRALPAK® AY-3R | 50 | 4.6 | 3 | Analytical |