



96 and 384 well Microplates Polypropylene

- Virgin polypropylene resin
- Biologically inert
- High chemical and thermal resistance (-196 °C to +121 °C)
- Manufactured in molds engineered specifically for polypropylene
- Conforms to ANSI/SBS - 1 (2004) Standards
- Footprint Dimensions to meet requirements for compatibility with equipment and robotics
- 96 well: Chimney-well design with flat, round and conical well bottoms having 390, 355 and 340 μ L maximum well volumes, respectively
- 384 well: Square round well top design with flat or conical well bottoms having 152 and 145 μ L maximum well volumes, respectively
- Sterile and pigmented versions

Free of detectable
DNase, RNase,
human DNA
non-pyrogenic

non-
cytotoxic

 All Greiner Bio-One plates can be barcoded to your specifications

GBO Cat. No.	Well Density	Bottom Style	Chimney Style	Color	Sterile	With Lid	Pack Qty	Case Qty
655201	96 well	Flat	Yes	Clear	No	No	10	100
655203	96 well	Flat	Yes	Red	No	No	10	100
655204	96 well	Flat	Yes	Blue	No	No	10	100
655205	96 well	Flat	Yes	Green	No	No	10	100
655206	96 well	Flat	Yes	Yellow	No	No	10	100
655209	96 well	Flat	Yes	Black	No	No	10	100
655261	96 well	Flat	Yes	Clear	Yes	No	10	100
650201	96 well	Round	Yes	Clear	No	No	10	100
650203	96 well	Round	Yes	Red	No	No	10	100
650204	96 well	Round	Yes	Blue	No	No	10	100
650205	96 well	Round	Yes	Green	No	No	10	100
650206	96 well	Round	Yes	Yellow	No	No	10	100
650207	96 well	Round	Yes	White	No	No	10	100
650209	96 well	Round	Yes	Black	No	No	10	100
650261	96 well	Round	Yes	Clear	Yes	No	10	100
651201	96 well	Conical	Yes	Clear	No	No	10	100
651203	96 well	Conical	Yes	Red	No	No	10	100
651204	96 well	Conical	Yes	Blue	No	No	10	100
651205	96 well	Conical	Yes	Green	No	No	10	100
651206	96 well	Conical	Yes	Yellow	No	No	10	100
651209	96 well	Conical	Yes	Black	No	No	10	100
651261	96 well	Conical	Yes	Clear	Yes	No	10	100
781201	384 well	Flat	No	Clear	No	No	10	100
781209	384 well	Flat	No	Black	No	No	10	100
781207	384 well	Flat	No	White	No	No	10	100
781280	384 well	Conical	No	Clear	No	No	10	100
781281	384 well	Conical	No	Clear	Yes	No	10	100