

## Injection of Liquid Samples

### Syringes For Manual Injection of Liquid Samples


Volume	5 $\mu\text{L}$	10 $\mu\text{L}$	25 $\mu\text{L}$	50 $\mu\text{L}$	100 $\mu\text{L}$	250 $\mu\text{L}$	500 $\mu\text{L}$
Model	75	701	702	705	710	725	750
Gauge	26s	26s	22s	22s	22s	22s	22
Cemented Needle (N) Syringe (point style 2)	87900	80300 80366 (6 pk)	80400	80500	80600	80700	80800
Removable Needle (RN) Syringe (point style 2)	87930	80330 80336 (6 pk)	80430	80530	80630	80730	80830
Replacement Needle 6 pk (point style 2)	7758-02	7758-02	7758-03	7758-03	7758-03	7779-03	7779-01

Note: Needles are 51 mm

### Zero Dead Volume Syringes

Volume	0.5 $\mu\text{L}$	1.0 $\mu\text{L}$	2.0 $\mu\text{L}$	5.0 $\mu\text{L}$
Model	7000.5	7001	7002	7105
Gauge	25	25	25	24
Knurled Hub (KH) Syringe (point style 2)	86259	80135	88411	88011
Standard Injection Spacer (25 mm)	86201	86201	86201	86201

Note: Needles are 70 mm



#### Tech Tip

Choose a syringe with a removable needle if working with halogenated solvents.