

Ordering information – LiChrospher® 100 RP-8 and RP-8 endcapped,
stainless steel cartridges LiChroCART®

Product	Ordering No.	Particle size	Dimension length	Dimension i.d.	Contents of one package
LiChrospher® 100 RP-8	1.50190.7087	5 µm	250 mm	2 mm	1 piece
LiChrospher® 100 RP-8	1.50177.7087	5 µm	250 mm	3 mm	1 piece
LiChrospher® 100 RP-8	1.50956.0001	5 µm	4 mm	4 mm	10 pieces
LiChrospher® 100 RP-8	1.50930.0001	5 µm	25 mm	4 mm	3 pieces
LiChrospher® 100 RP-8	1.50228.7087	5 µm	55 mm	4 mm	1 piece
LiChrospher® 100 RP-8	1.50986.0001	5 µm	75 mm	4 mm	3 pieces
LiChrospher® 100 RP-8	1.50822.0001	5 µm	125 mm	4 mm	1 piece
LiChrospher® 100 RP-8	1.50942.0001	5 µm	125 mm	4 mm	3 pieces
LiChrospher® 100 RP-8	1.50832.0001	5 µm	250 mm	4 mm	1 piece
LiChrospher® 100 RP-8	1.50982.0001	5 µm	250 mm	4 mm	3 pieces
LiChrospher® 100 RP-8	1.50634.0001	5 µm	100 mm	4.6 mm	1 piece
LiChrospher® 100 RP-8	1.50635.0001	5 µm	150 mm	4.6 mm	1 piece
LiChrospher® 100 RP-8	1.50636.0001	5 µm	250 mm	4.6 mm	1 piece
LiChrospher® 100 RP-8	1.51443.7087	5 µm	125 mm	10 mm	1 piece
LiChrospher® 100 RP-8	1.50842.0001	10 µm	250 mm	4 mm	1 piece
LiChrospher® 100 RP-8	1.51448.7088	10 µm	100 mm	4.6 mm	1 piece
LiChrospher® 100 RP-8	1.50179.7088	10 µm	250 mm	10 mm	1 piece
LiChrospher® 100 RP-8 endcapped	1.50175.7092	5 µm	125 mm	3 mm	1 piece
LiChrospher® 100 RP-8 endcapped	1.50177.7092	5 µm	250 mm	3 mm	1 piece
LiChrospher® 100 RP-8 endcapped	1.50961.0001	5 µm	4 mm	4 mm	10 pieces
LiChrospher® 100 RP-8 endcapped	1.50827.0001	5 µm	125 mm	4 mm	1 piece
LiChrospher® 100 RP-8 endcapped	1.50837.0001	5 µm	250 mm	4 mm	1 piece
LiChrospher® 100 RP-8 endcapped	1.50637.0001	5 µm	100 mm	4.6 mm	1 piece
LiChrospher® 100 RP-8 endcapped	1.51442.7092	5 µm	125 mm	4.6 mm	1 piece
LiChrospher® 100 RP-8 endcapped	1.50638.0001	5 µm	150 mm	4.6 mm	1 piece
LiChrospher® 100 RP-8 endcapped	1.50639.0001	5 µm	250 mm	4.6 mm	1 piece
LiChrospher® 100 RP-8 endcapped	1.50234.7091	10 µm	55 mm	2 mm	1 piece
LiChrospher® 100 RP-8 endcapped	1.50236.7091	10 µm	55 mm	3 mm	1 piece
LiChrospher® 100 RP-8 endcapped	1.50847.0001	10 µm	250 mm	4 mm	1 piece
LiChrospher® 100 RP-8 endcapped	1.51432.7091	10 µm	150 mm	4.6 mm	1 piece
LiChrospher® 100 RP-8 endcapped	1.51431.7091	10 µm	250 mm	4.6 mm	1 piece
LiChrospher® 100 RP-8 endcapped	1.50179.7091	10 µm	250 mm	10 mm	1 piece

The LiChroCART® columns (75, 100, 125, 150 and 250 mm length) in the list above (4 or 4.6 mm i.d.) require part number 1.51486.0001 manu-CART® cartridge column holder, which can be used to hold one cartridge column with or without a 4–4 mm guard column. LiChroCART® columns 250–10 mm require part number 1.51419.0001 manu-CART® 10. Additional dimensions available as customized packings see page 366.