

## Ordering information – Superspher®, stainless steel cartridges LiChroCART®

Product	Ordering No.	Particle size	Dimension length	Dimension i.d.	Contents of one package
Superspher® 100 RP-18	1.50204.0001	4 µm	10 mm	2 mm	3 pieces
Superspher® 100 RP-18	1.50200.0001	4 µm	125 mm	2 mm	1 piece
Superspher® 100 RP-18 (3 HPLC columns from 3 batches)	1.50200.1003	4 µm	125 mm	2 mm	3 pieces
Superspher® 100 RP-18	1.50190.7137	4 µm	250 mm	2 mm	1 piece
Superspher® 100 RP-18	1.50792.0001	4 µm	125 mm	3 mm	1 piece
Superspher® 100 RP-18 (3 HPLC columns from 3 batches)	1.50792.1003	4 µm	125 mm	3 mm	3 pieces
Superspher® 100 RP-18	1.51299.0001	4 µm	250 mm	3 mm	1 piece
Superspher® 100 RP-18 (3 HPLC columns from 3 batches)	1.51299.1003	4 µm	250 mm	3 mm	3 pieces
Superspher® 100 RP-18	1.16039.0001	4 µm	25 mm	4 mm	3 pieces
Superspher® 100 RP-18	1.50980.0001	4 µm	75 mm	4 mm	3 pieces
Superspher® 100 RP-18	1.16051.0001	4 µm	125 mm	4 mm	1 piece
Superspher® 100 RP-18 (3 HPLC columns from 3 batches)	1.16051.1003	4 µm	125 mm	4 mm	3 pieces
Superspher® 100 RP-18	1.16056.0001	4 µm	250 mm	4 mm	1 piece
Superspher® 100 RP-18 (3 HPLC columns from 3 batches)	1.16056.1003	4 µm	250 mm	4 mm	3 pieces
Superspher® 100 RP-18 endcapped	1.50198.0001	4 µm	125 mm	2 mm	1 piece
Superspher® 100 RP-18 endcapped	1.50193.0001	4 µm	250 mm	2 mm	1 piece
Superspher® 100 RP-18 endcapped	1.51909.0001	4 µm	125 mm	3 mm	1 piece
Superspher® 100 RP-18 endcapped	1.51910.0001	4 µm	250 mm	3 mm	1 piece
Superspher® 100 RP-18 endcapped	1.16869.0001	4 µm	25 mm	4 mm	3 pieces
Superspher® 100 RP-18 endcapped	1.50302.7138	4 µm	30 mm	4 mm	1 piece
Superspher® 100 RP-18 endcapped	1.16855.0001	4 µm	125 mm	4 mm	1 piece
Superspher® 100 RP-18 endcapped	1.16858.0001	4 µm	250 mm	4 mm	1 piece
Superspher® 100 RP-18 endcapped	1.51448.7138	4 µm	100 mm	4.6 mm	1 piece
Superspher® 100 RP-18 endcapped	1.51442.7138	4 µm	125 mm	4.6 mm	1 piece
Superspher® 100 RP-18 endcapped	1.51432.7138	4 µm	150 mm	4.6 mm	1 piece
Superspher® 100 RP-18 endcapped	1.51431.7138	4 µm	250 mm	4.6 mm	1 piece

The LiChroCART® columns (75, 125, 150 and 250 mm length) in the list on the left (2, 3, 4 and 4.6 mm i.d.) require part number 1.51486.0001 manu-CART® cartridge column holder, which can be used to hold one cartridge column with or without a 4–4 mm guard column. Additional dimensions available as customized packings see page 366. As guard column we recommend LiChroCART® 4–4 LiChrospher® guard cartridges.