

Stainless Steel Forceps

Make filter handling easy



- ▶ Tips have a flat, smooth surface to avoid membrane filter damage.
- ▶ Polypropylene finger grips provide a comfortable and secure hold.
- ▶ Choose traditional black or multi-colored finger grips. Bright colors make forceps easy to identify, track, and see on the lab bench.

Applications

- ▶ Ideal for handling and moving membrane to and from filter holders and Petri dishes.

Specifications

Materials of Construction

Stainless steel with polypropylene finger grips

Sterilization

Provided non-sterile. "Flame" the tips prior to use, or autoclave at 121 - 123 °C (250 - 253 °F) at 1.0 bar (100 kPa, 15 psi) for 15 - 20 min. Do not autoclave in aluminum foil.

Ordering Information

Stainless Steel Forceps

Part Number	Description	Pkg
51147	Stainless steel forceps, black grips	1/pkg
4690	Stainless steel forceps, multi-colored grips (1 each of orange, blue, chartreuse)	3/pkg

Petri Dishes

Unique patented design for easy handling and storage



- ▶ Opens easily with one hand, yet closes to a tight seal.
- ▶ Uses less space on the lab bench or in the incubator with easy stacking base.
- ▶ Gamma-irradiated. No EtO residuals to impede microbial growth.
- ▶ Available with or without absorbent pads.
- ▶ Bulk packaging offers additional value.

Applications

- ▶ Ideal for microbiological analysis when performing the Membrane Filter (MF) Technique.
- ▶ Petri dishes with absorbent pads can be used with broth media, or users can pour agar into dishes without absorbent pads.

Specifications

Petri Dish

Material of Construction

Polystyrene

Dimensions

Height: 9.0 mm (0.35 in.)
Diameter: 50.0 mm (1.97 in.)

Filter Size

Accepts 47 mm membrane filter

Gamma-Irradiated

Validated dose 15 - 30 kGy

Absorbent Pads

Composition

Cellulose

Typical Thickness

0.9 mm (35 mils)

Ordering Information

Petri Dishes

Part Number	Description	Pkg
7242	Petri dishes, without absorbent pads	100/pkg
7232	Petri dishes, bulk pack, without absorbent pads	500/pkg
7245	Petri dishes, with absorbent pads	100/pkg