

## Kalrez O-Ring Kits for Sampling Accessories

| Description  | Part No.        |
|--|-----------------|
| Nebulizer, End Cap, and Burner Head O-Ring Kit for Plastic Nebulizer for PinAAcle AA.<br>Includes: <b>09200272</b> (1), <b>09200275</b> (4), <b>09200276</b> (2), <b>09200277</b> (1), <b>09200278</b> (1), <b>09200369</b> (2)                              | <b>N3160157</b> |
| Nebulizer, End Cap, and Burner Head O-Ring Kit for Stainless Steel Nebulizer for PinAAcle AA.<br>Includes: <b>09200271</b> (1), <b>09200273</b> (2), <b>09200274</b> (1), <b>09200275</b> (4), <b>09200276</b> (2), <b>09200277</b> (1), <b>09200278</b> (1) | <b>N3160156</b> |
| Nebulizer O-Ring Kit for Stainless Steel Nebulizer for PinAAcle AA. Includes: <b>09200271</b> (1), <b>09200273</b> (2), <b>09200274</b> (1), <b>09200276</b> (2)   | <b>N3160192</b> |
| Nebulizer O-Ring Kit, fits Stainless Steel Nebulizers <b>03030358</b> and <b>03030352</b> . Includes: <b>09200271</b> (1), <b>09200273</b> (2), <b>09200274</b> (1)  | <b>N2010247</b> |
| Nebulizer O-Ring Kit for Plastic Nebulizer for PinAAcle AA. Includes: <b>09200272</b> (1), <b>09200276</b> (2), <b>09200369</b> (2)  | <b>N3160193</b> |
| Nebulizer O-Ring Kit for AA200/400 Nebulizer <b>N2010250</b> . Includes: <b>09200272</b> (1), <b>09200276</b> (2), <b>09200369</b> (2)   | <b>N2010246</b> |
| End Cap O-Ring Kit for PinAAcle AA. Includes: <b>09200275</b> (4), <b>09200277</b> (1)   | <b>N3160194</b> |
| End Cap O-Ring Kit for AAnalyst 700/800 End Cap, <b>00570984</b> . Includes: <b>09200277</b> (1), <b>09200341</b> (2), <b>09200342</b> (2)   | <b>N2010245</b> |

## PinAAcle Nebulizers and Assemblies Quick Reference

| Description  | Type                | Needle           | O-Ring Type | Complete Assembly     | End Cap Assembly | Replacement Needle Assembly |
|--|---------------------|------------------|-------------|-----------------------|------------------|-----------------------------|
| Plastic High Sensitivity (standard)  | Aqueous             | Platinum/Iridium | Viton       | <b>N3160144</b>       | <b>N3160511</b>  | <b>B3150525</b>             |
| Plastic High Sensitivity   | Aqueous or Organics | Platinum/Iridium | Kalrez      | <b>N3160159</b>       | <b>N3160189</b>  | <b>N3160154</b>             |
| Plastic High Sensitivity (suitable for aqua regia in addition to organics) | Aqueous or Organics | Tantalum         | Kalrez      | <b>N3160152</b>       | <b>N3160189</b>  |                             |
| Stainless Steel (standard)   | Aqueous             | Stainless        | Viton       | <b>N3160143</b>       | <b>N3160505</b>  | <b>03030354</b>             |
| Stainless Steel  | Aqueous or Organics | Stainless        | Kalrez      | <b>N3160158</b>       | <b>N3160188</b>  | <b>N3160153</b>             |
| Stainless Steel  | Aqueous             | Plastic          | Viton       |                       | <b>N3160505</b>  | <b>03030402</b>             |
| Stainless Steel  | Aqueous             | Platinum Alloy   | Viton       |                       | <b>N3160505</b>  | <b>03030133</b>             |
| Description  | Type                |                  |             | Burner Drain Assembly | Drain Bottle     |                             |
| Drain Assembly   | Aqueous             |                  |             | <b>N0400058</b>       | <b>09200486</b>  |                             |
| Drain Assembly   | Organics            |                  |             | <b>N3160683</b>       | <b>09200486</b>  |                             |

## AAAnalyst Nebulizers and Assemblies Quick Reference

| Description                    | Type                | Needle                 | O-Ring Type | Complete Assembly | End Cap Assembly | Replacement Needle Assembly |
|--------------------------------|---------------------|------------------------|-------------|-------------------|------------------|-----------------------------|
| High-Sensitivity Plastic Body  | Aqueous or Organics | Tantalum               | Kalrez      | <b>N2010250</b>   |                  | <b>N3160146</b>             |
| Plastic Body                   | Aqueous             | Platinum/Iridium       | Viton       | <b>N3150188</b>   | <b>N0370392</b>  | <b>N3150525</b>             |
| High-Sensitivity Plastic Body  | Aqueous or Organic  | Platinum/Iridium       | Kalrez      | <b>B0505480</b>   | <b>N0370392</b>  | <b>B0501022</b>             |
| High-Sensitivity Plastic Body  | Aqueous or Organic  | Tantalum               | Kalrez      | <b>B0505590</b>   | <b>N0370392</b>  |                             |
| Stainless Steel Body           | Aqueous             | Stainless Steel        | Viton       | <b>03030352</b>   | <b>00570984</b>  | <b>03030354</b>             |
| Corrosion Resistant Metal Body | Aqueous or Organic  | Tantalum               | Kalrez      | <b>03030404</b>   | <b>00570984</b>  | <b>03030402</b>             |
| Platinum Alloy Nebulizer       | Aqueous             | Platinum/Iridium       | Kalrez      | <b>03030299</b>   | <b>00570984</b>  | <b>03030433</b>             |
| Universal GemTip               | Aqueous             | Plastic w/Sapphire Tip | Viton       | <b>N0370394</b>   | <b>N0370392</b>  | <b>N0371520</b>             |
| High-Sensitivity GemTip        | Aqueous or Organic  | Plastic w/Sapphire Tip | Kalrez      | <b>N0370393</b>   | <b>N0370392</b>  | <b>N0371522</b>             |

## Nebulizer Accessories

| Description                     | Part No.        |
|---------------------------------|-----------------|
| Capillary Tubing                | <b>09908265</b> |
| Nebulizer Cleaning Wires Pkg. 5 | <b>03030135</b> |
| Nebulizer Grease                | <b>03030405</b> |



Nebulizer Capillary Tubing



Nebulizer Cleaning Wires



Nebulizer Grease