


# Packed Traps for PerkinElmer Thermal Desorbers

## Air Monitoring Trap

Low flow trap packed with carbonaceous sorbents suitable for ozone precursor and air toxin monitoring.



## TurboMatrix Thermal Desorption Cold Trap Supplies

Trap supplies from PerkinElmer, the market leader in thermal desorption, will provide exceptional analytical performance. Used for U.S. EPA Method TO17, the PerkinElmer standard trap, packed with Tenax™, on the TurboMatrix™ Thermal Desorber, will improve productivity and trapping capacity. The TurboMatrix air monitoring trap is packed with carbonaceous sorbents suitable for ozone precursor and air toxics monitoring.

Description	Qty.	Part No.
<b>Cold Traps for TurboMatrix</b>		
Air Monitoring Trap	1	<b>M0413628</b>
Empty Trap	1	<b>M0413627</b>
Tenax™ TA 60/80 Packed Trap	1	<b>M0413535</b>
Carbopack® C Packed Trap	1	<b>N6200150</b>
<b>Cold Traps for ATD 400</b>		
Air Monitoring Trap	1	<b>L4275108</b>
Trap Nuts (2 required)	1	<b>L4275009</b>
Trap Tube Low Flow, Empty (Narrow-bore at one end. Allows minimum gas flow during trap desorption)	1	<b>L4275107</b>
Tenax™ TA 60/80 Packed Trap	1	<b>L4275089</b>
Empty Trap	1	<b>L4271106</b>



## Cold Traps Fittings and Accessories

Description	Qty.	Part No.
Graphite Ferrule	2	<b>L4271187</b>
SilTite™ Ferrule (GC/MS) 0.4 mm	10	<b>N9306093</b>
SilTite™ Ferrule (GC/MS) 0.5 mm	10	<b>N9306094</b>
SilTite™ Ferrule (GC/MS) 0.8 mm	10	<b>N9306095</b>
SilTite™ Nuts	5	<b>N9306096</b>
PTFE Ferrule	10	<b>L4275110</b>
Graphite/Vespel® Ferrule for ATD 400	5	<b>L1003027</b>
Valco® Graphite/Vespel® Ferrule for ATD 400	5	<b>L1003028</b>
Trap Filter Disk		<b>L1003030</b>
Trap Packing Disk	20	<b>L4271290</b>
Trap Packing Retaining Spring	5	<b>N6301054</b>
Quartz Wool (Untreated) 5 g		<b>N6102354</b>
Cold Trap Packing Tool		<b>L4271203</b>
Cold Trap Removal Tool		<b>L4271205</b>
Regulator 0–60 psig		<b>N6101474</b>
Backflush Nozzle for ATD 400		<b>L4275072</b>
Internal Standard Injection Accessory for ATD 400		<b>L4270010</b>
Liquid Nitrogen Accessory for ATD 400		<b>L4270009</b>
Gauze Loading Rig		<b>L4070023</b>
Replacement Plastic Plunger for Gauze Loading Rig		<b>L4071151</b>