



## Quick-Change Torch Module with Cyclonic Spray Chamber and Concentric Nebulizer

ICP Model	Part No.
<b>Optima 8x00</b> Includes:	<b>N0780607</b>
MEINHARD Type K1 Nebulizer 1 mL/Min	<b>N0777707</b>
1 Slot Quartz Torch	<b>N0780130</b>
Alumina Injector 2.0 mm ID	<b>N0781014</b>
Polyethylene Tubing, 10 ft.	<b>09908265</b>
Adjustable Torch Assembly	<b>N0770600</b>
Cyclonic Spray Chamber w/Adjustable Mount	<b>N0770605</b>
Tubing Assembly	<b>N0770610</b>
<b>Optima 2x00/4x00/5x00/7x00</b> Includes:	<b>N0770607</b>
MEINHARD Type C1 Nebulizer 1 mL/Min	<b>00472022</b>
Polyethylene Tubing, 10 ft.	<b>09908265</b>
1 Slot Quartz Torch	<b>N0770338</b>
Adjustable Torch Assembly	<b>N0770600</b>
Cyclonic Spray Chamber w/Adjustable Mount	<b>N0770605</b>
Tubing Assembly	<b>N0770610</b>
Alumina Injector 2.0 mm	<b>N0775177</b>

## Quick-Change Torch Module for HF Analysis

ICP Model	Part No.
<b>Optima 4300V/5300V/7300V</b> Includes:	<b>N0770911</b>
GemCone High Dissolved Solids Nebulizer	<b>N0690670</b>
Alumina Sample Injector 2.0 mm ID	<b>N0695362</b>
Quick Disconnect Tubing	<b>N0770336</b>
3 Slot Quartz Torch	<b>N0771500</b>
Tracey TFE Spray Chamber Includes: Fittings, O-Rings, and Spares	<b>N0777497</b>

## Quick-Change Torch Module for Oil Analysis

ICP Model	Part No.
<b>Optima 4300V/5300V/7300V</b> Includes:	<b>N0770910</b>
GemCone Low-Flow Nebulizer	<b>N0770358</b>
Alumina Injector 1.2 mm	<b>N0771531</b>
Quick Disconnect Tubing	<b>N0770336</b>
3 Slot Quartz Torch	<b>N0771500</b>
Baffled Cyclonic Spray Chamber Includes: Fittings, O-Rings, and Spares	<b>N0771501</b>



## Quick-Change Torch Module with Scott Spray Chamber and GemTip Cross-Flow II Nebulizer

ICP Model	Part No.
<b>Optima 8x00</b> Includes:	<b>N0780606</b>
1 Slot Quartz Torch	<b>N0780130</b>
GemTip Cross-Flow II Nebulizer	<b>N0780546</b>
Alumina Injector 2.0 mm ID	<b>N0781014</b>
Scott Spray Chamber	<b>N0770357</b>
<b>Optima 2x00/4x00/5x00/7x00 DV</b> Includes:	<b>N0770606</b>
1 Slot Quartz Torch	<b>N0770338</b>
GemTip CrossFlow II Nebulizer	<b>N0770546</b>
Adjustable Torch Assembly	<b>N0770600</b>
Alumina Injector 2.0 mm	<b>N0775177</b>
Scott Spray Chamber	<b>N0770357</b>

## Mounting Block

Mounting block only. Does not include torch, injector, spray chamber or nebulizer.



ICP Model	Part No.
Optima 2x00/4x00/5x00/7x00/8x00	<b>N0770600</b>

## Quick-Change Torch Module Support Stand

A free-standing support stand to hold your spare Quick-Change Torch Module.



ICP Model	Part No.
Optima 2x00/4x00/5x00/7x00 DV/8x00	<b>N0770609</b>

## End Cap for Concentric Glass, GemCone and Similar Nebulizers

End cap to fit any concentric nebulizer onto a Scott Spray Chamber.



ICP Model	Part No.
Optima 2x00/3x00/4x00/5x00/7x00/8x00 / Avio 200/500	<b>N0680504</b>