



Part No. **N0771500**



Part No. **N0772006**

## Organics Kits

ICP Model	Qty.	Part No.
<b>Avio 200/500</b> Includes:		<b>N0790610</b>
3 Slot Quartz Torch	2	<b>N0790249</b>
1.2 mm ID Alumina Injector (1.2 mm entire length of tube)	2	<b>N0791182</b>
Baffled Cyclonic Spray Chamber	1	<b>N0791352</b>
GemCone™ Low-Flow Nebulizer	1	<b>N0770358</b>
Cyclonic Spray Chamber Adapter	1	<b>N0790141</b>
Nebulizer Spacer	1	<b>N0791438</b>
<b>Optima 8x00</b> Includes:		<b>N0780608</b>
3 Slot Quartz Torch	2	<b>N0780132</b>
1.2 mm ID Alumina Injector	2	<b>N0781007</b>
Baffled Cyclonic Spray Chamber 4 mm	1	<b>N0776090</b>
GemCone™ Low-Flow Nebulizer	1	<b>N0770358</b>
Cyclonic Spray Chamber Adapter	1	<b>N0770614</b>
Injector Adapter Washer	1	<b>N0776014</b>
Injector Adapter	1	<b>N0776091</b>
Injector Adapter Nut	1	<b>N0776025</b>
O-Ring	1	<b>09902207</b>
O-Ring 8.74 ID 1.78 WD	2	<b>09921028</b>
<b>Optima 2x00/4x00/5x00/7x00 DV</b> Includes:		<b>N0770420</b>
3 Slot Quartz Torch	1	<b>N0772005</b>
1.2 mm ID Alumina Injector	1	<b>N0776093</b>
Baffled Cyclonic Spray Chamber 4 mm	1	<b>N0776090</b>
GemCone™ Low-Flow Nebulizer	1	<b>N0770358</b>
Cyclonic Spray Chamber Adapter	1	<b>N0770614</b>
Injector Adapter Washer	1	<b>N0776014</b>
Injector Adapter Nut	1	<b>N0776025</b>
Injector Adapter	1	<b>N0776091</b>
O-Ring	1	<b>09902207</b>
O-Ring 8.74 ID 1.78 WD	2	<b>09921028</b>
<b>Optima 4300 V/5300 V/7300 V</b> Includes:		<b>N0770732</b>
3 Slot Quartz Torch	1	<b>N0771500</b>
Quartz Bonnet	1	<b>N0775289</b>
1.2 mm ID Alumina Injector	1	<b>N0771531</b>
Baffled Cyclonic Spray Chamber	1	<b>N0771501</b>
GemCone™ Low-Flow Nebulizer	1	<b>N0770358</b>
Quartz Torch Copper Ignitor Tape	120	<b>N0775297</b>

## Oxygen Addition Kits

ICP Model	Qty.	Part No.
<b>Oxygen Addition Kit</b> <b>Avio 200/500</b> Includes:		<b>N0790412</b>
Baffled Cyclonic Spray Chamber 4 mm for Oils with Oxygen Port	1	<b>N0770389</b>
Tygon Tubing	25 ft.	<b>02506516</b>
Insert	4	<b>09903094</b>
Tube Fitting	2	<b>09903150</b>
Mass Flow Controller	1	<b>N0772006</b>
<b>ICP Model</b>	<b>Qty.</b>	<b>Part No.</b>
<b>DV Oxygen Addition Kit</b> <b>Optima 2x00/4x00/5x00/7x00/8x00 DV</b> Includes:		<b>N0770689</b>
Baffled Cyclonic Spray Chamber 4 mm for Oils with Oxygen Port	1	<b>N0777731</b>
Tygon Tubing	25 ft.	<b>02506516</b>
Insert	4	<b>09903094</b>
Tube Fitting	2	<b>09903150</b>
Mass Flow Controller	1	<b>N0772006</b>
Adapter for Cyclonic Spray Chamber	1	<b>N0770614</b>
<b>ICP Model</b>	<b>Qty.</b>	<b>Part No.</b>
<b>V Oxygen Addition Kit</b> <b>Optima 4300 V/5300 V/7300 V</b> Includes:		<b>N0770965</b>
Torch Adapter with Auxiliary Oxygen Port	1	<b>N0774070</b>
Tygon Tubing Per l/0.3 m	25 ft.	<b>02506516</b>
Insert	4	<b>09903094</b>
Tube Fitting	2	<b>09903150</b>
Mass Flow Controller	1	<b>N0772006</b>