



## Carbon Fiber Supported Autosampler Probes for ST Nebulizers


Self-aspiration rate @ 1 L/m Ar

Description	Part No.
0.25 mm ID, Capillary Length 80 cm, Self-Aspiration Rate (@ 1 L/m Ar) 100 µL/min	<b>N0777223</b>
0.50 mm ID, Capillary Length 80 cm, Self-Aspiration Rate (@ 1 L/m Ar) 400 µL/min	<b>N0777225</b>
0.25 mm ID, Capillary Length 150 cm, Self-Aspiration Rate (@ 1 L/m Ar) 100 µL/min	<b>N0777826</b>
0.50 mm ID, Capillary Length 150 cm, Self-Aspiration Rate (@ 1 L/m Ar) 100 µL/min	<b>N0777827</b>
Reset Probe Weight Kit. PTFE coated Stainless Steel Weight with Black PEEK Retainer Assembly	<b>N0786043</b>

## PerkinElmer S10 Autosampler Probes for ST Nebulizers

Description	Part No.
PTFE Flame Sampling Probe Assembly with Carbon Fiber Support and Female Fitting, 0.25 mm ID for ICP-MS. Requires but does not include <b>N8145017</b>	<b>N8145139</b>
PTFE Flame Sampling Probe Assembly with Carbon Fiber Support and Female Fitting, 0.5 mm ID for ICP/ICP-MS. Requires but does not include <b>N8145017</b>	<b>N8145140</b>
PTFE Flame Sampling Probe Assembly with Carbon Fiber Support and Male Screw Fitting, 0.55 mm ID for AA/ICP/ICPMS. For ICP/ICPMS requires but does not include <b>N8145018*</b> . For AA with a stainless steel nebulizer, the male connect should be removed so that tubing can be directly connected to the nebulizer, no additional tubing is required	<b>N2024132</b>
PTFE Flame Sampling Probe Assembly and Male Screw Fitting, 0.8 mm ID for AA. Requires but does not include <b>B0172053*</b>	<b>N0777434</b>
PTFE Flame Sampling Probe Assembly with Carbon Fiber Support and Male Screw Fitting, 1.0 mm ID for FIAS*	<b>N2024131</b>

\*Keep spare ferrules on hand, **N8145237** Qty. 10



## PFA MicroFlow Nebulizers with Integrated Capillaries for NexION 300/350

Measured self-aspiration rate @ 1 L/m Ar

Description	Uptake Rate	Part No.
PFA-20 PFA Nebulizer	20 µL/min	<b>N8122350</b>
PFA-50 PFA Nebulizer	50 µL/min	<b>N8145102</b>
PFA-100 PFA Nebulizer	100 µL/min	<b>N8122351</b>
PFA-200 PFA Nebulizer	200 µL/min	<b>N8145103</b>
PFA-400 PFA Nebulizer	400 µL/min	<b>N8145104</b>



## PFA MicroFlow Nebulizers with Integrated Autosampler Probes for NexION 300/350

Measured self-aspiration rate @ 1 L/m Ar

Description	Uptake Rate	Part No.
PFA-20 PFA (Carbon Fiber Support)	20 µL/min, length 80 cm	<b>N0777214</b>
PFA-50 PFA (Carbon Fiber Support)	50 µL/min, length 80 cm	<b>N0777215</b>
PFA-100 PFA (Carbon Fiber Support)	100 µL/min, length 100 cm	<b>N0777216</b>
PFA-200 PFA (Carbon Fiber Support)	200 µL/min, length 100 cm	<b>N0777217</b>
PFA-400 PFA (Carbon Fiber Support)	400 µL/min, length 100 cm	<b>N0777218</b>
PFA-20 PFA (Ultem Support)	20 µL/min, length 100 cm	<b>N8145141</b>
PFA-50 PFA (Ultem Support)	50 µL/min, length 100 cm	<b>N8145142</b>
PFA-100 PFA (Ultem Support)	100 µL/min, length 100 cm	<b>N8145143</b>
PFA-200 PFA (Ultem Support)	200 µL/min, length 100 cm	<b>N8145144</b>
PFA-400 PFA (Ultem Support)	400 µL/min, length 100 cm	<b>N8145275</b>

## PFA Nebulizers with Integrated Capillaries for NexION 1000/2000

Measured self-aspiration rate @ 1 L/m Ar

Description	Uptake Rate	Part No.
PFA-20 PFA Nebulizer	20 µL/min	<b>N8152582</b>
PFA-50 PFA Nebulizer	50 µL/min	<b>N8152583</b>
PFA-100 PFA Nebulizer	100 µL/min	<b>N8152584</b>
PFA-200 PFA Nebulizer	200 µL/min	<b>N8152585</b>
PFA-400 PFA Nebulizer	400 µL/min	<b>N8152586</b>

## PFA Nebulizers with Integrated Autosampler Probes for NexION 1000/2000

Measured self-aspiration rate @ 1 L/m Ar

Description	Uptake Rate	Part No.
PFA-20 PFA (Carbon Fiber Support)	20 µL/min, length 80 cm	<b>N8152587</b>
PFA-50 PFA (Carbon Fiber Support)	50 µL/min, length 80 cm	<b>N8152588</b>
PFA-100 PFA (Carbon Fiber Support)	100 µL/min, length 100 cm	<b>N8152589</b>
PFA-200 PFA (Carbon Fiber Support)	200 µL/min, length 100 cm	<b>N8152590</b>
PFA-400 PFA (Carbon Fiber Support)	400 µL/min, length 100 cm	<b>N8152591</b>
PFA-20 PFA (Ultem Support)	20 µL/min, length 100 cm	<b>N8152592</b>
PFA-50 PFA (Ultem Support)	50 µL/min, length 100 cm	<b>N8152593</b>
PFA-100 PFA (Ultem Support)	100 µL/min, length 100 cm	<b>N8152594</b>
PFA-200 PFA (Ultem Support)	200 µL/min, length 100 cm	<b>N8152420</b>
PFA-400 PFA (Ultem Support)	400 µL/min, length 100 cm	<b>N8152595</b>