

FID Miscellaneous Accessories*

| Description | Part No. |
|---|----------|
| Auto-Ignite Nozzle Replacement Tool for removing nozzle from auto-ignite FID body | N6103188 |
| Auto-Ignite FID Catalytic Reactor Accessory* | N6120161 |
| Capillary Column Adapter for capillary column use with the FID 1/8 to 1/16 in. detector adapter | N6120020 |
| Ceramic Column Cutter | N9301376 |
| Cotton Applicators | N9301272 |
| Detector Cover (Brown Color) | N6103225 |
| Detector Cover (White Color) | N6103151 |
| Eraser Brush-Pencil | 09923078 |
| Jet Assembly* | N6100194 |
| Jet Replacement Tool 1/4 in. nut driver for removing jet from auto-ignite FID body | N6101297 |
| Quartz Wool | N6102354 |
| Replacement Stainless Steel Glow Plug for Auto-ignite FID | N6103089 |
| Septa, low bleed (pkg. 50) | N9302972 |
| 1/4 in. Packed Column Adapter for use with 1/4 in. packed columns. 1/8 to 1/4 in. adapter fits both injector and detector ends | 00080100 |

*Not compatible with the FID detector on the Clarus 690 and 590 GCs
Refer to the Clarus 590/690 consumables reference guide for more details

Catalytic Reactor Accessory

Auto-Ignite FID with PPC Add-On Kit

The catalytic reactor converts CO and CO₂ to methane conveniently and efficiently. The lower detection limit is extended to well below 0.1 ppm.

The catalytic reactor consists of a special catalytic reactor base which replaces the FID base. The reactor body contains a quantity of catalyst held in place by a quartz wool plug.

Kits include base assembly, reactor tube, and instructions.

External Igniter FID Catalytic Reactor Accessory

| Description | Part No. |
|---|----------|
| External Igniter FID Catalytic Reactor Accessory* | N6120070 |

Auto-Ignite FID Catalytic Reactor Accessory


| Description | Part No. |
|--|----------|
| Auto-Ignite FID Catalytic Reactor Accessory* | N6120161 |

*Not compatible with the FID detector on the Clarus 690 and 590 GCs
Refer to the Clarus 590/690 consumables reference guide for more details

Catalytic Reactor Replacement Parts

| Description | Part No. |
|--------------|----------|
| Catalyst | N9302698 |
| Jet Assembly | N6100194 |
| Quartz Wool* | N6102354 |

* Product available in the USA only



Miscellaneous Accessories

External Ignite Nozzle Replacement Tool

For removing nozzle from External Ignite FID body.

| Description | Part No. |
|---|----------|
| External Ignite Nozzle Replacement Tool | N6103188 |

Hydrogen Regulator Replacement Kit

| Description | Part No. |
|------------------------------------|----------|
| Hydrogen Regulator Replacement Kit | N6100289 |

Hydrogen/Air Replacement Needle Valve

| Description | Part No. |
|---------------------------------------|----------|
| Hydrogen/Air Replacement Needle Valve | N6101412 |

Igniter Assembly

Complete replacement igniter assembly with glow plug.

| Description | Part No. |
|------------------|----------|
| Igniter Assembly | N6100016 |

Replacement Glow Plug for External Ignite FID

| Description | Part No. |
|---|----------|
| Replacement Glow Plug for External Ignite FID | 00091279 |

FID Flow Measurement Adapter

| Description | Part No. |
|------------------------------|----------|
| FID Flow Measurement Adapter | N6101345 |