

UPS Systems

On-Line Conditioned Uninterruptible Power Supply Systems

With the Security Plus Series, you get much more protection and a higher comfort level than you get with most other UPS systems. The Security Plus Series also provides complete power conditioning and, because the Security Plus Series features on-line inverter design, added peace of mind. And regardless of input fluctuations, the Security Plus Series insures that the output remains continuous and regulated.




Instrument Model	Description	Frequency (Hz)	Part No.
Atomic Absorption			
PinAAcle 500	1440 VA	60	N9306750
PinAAcle 500	1440 VA	50	N0777798
PinAAcle 900F/AAAnalyst 100/200/300/400	800 VA	60	N0777681
PinAAcle 900F/AAAnalyst 100/200/300/400	800 VA	50	N0777689
AAAnalyst 700/HGA 900	5.2 kVA	50/60	N9308150
PinAAcle 900 T/H/Z	10 kVA	50/60	N9306757
AAAnalyst 600/800	10 kVA	50/60	N0777613
Optima™ 2x00/4x00/5x00/7x00/8x00	3.8 kVA	60	N9307512
Optima 2x00/4x00/5x00/7x00/8x00	3.6 kVA	50	N9307522
Inductively Coupled Plasma-OES			
Optima 2x00/4x00/5x00/7x00/8x00	5.2 kVA	50/60	N0777511
Inductively Coupled Plasma-MS			
NexION 300/350/1000/2000	8.0 kVA	50/60	N0777682
ELAN	12 kVA	50/60	N0777719

UPS Systems include a POWERVAR three-year warranty (two-year on batteries) on parts and labor.

Extended Run Battery Cabinets

Description	Shipping Weight	Part No.
Extended Run Battery Cabinet 30 min additional run time battery for Optima/Avio. 208/240 Volt, 8 Amp input for internal charger Compatible with: N0777511 or N9308150	413 lbs. (187 kg)	N0777727
Extended Run Battery Cabinet 60 min additional run time battery for Optima/Avio. 208/240 Volt, 8 Amp input for internal charger Compatible with: N0777511 or N9308150	497 lbs. (225 kg)	N0777728
Extended Battery Cabinet 45 min additional run time battery for NexION. 208/240 Volt input for internal charger Compatible with: N0777682	746 lbs. (338 kg)	N0777729



The ABC's of power conditioning