

Injector Assemblies and Support Adapters



Part No. **N8122358**

PFA Injector Assemblies

Description	Recommended Use	Part No.
PFA/Quartz Torch Injector Assembly		
PFA/Quartz Twist Type	1.5 mm ID	N8122394
PFA/Quartz Cassette	1.5 mm ID	N8122413
PFA/Sapphire Torch Injector Assembly		
PFA/Sapphire Twist Type	1.8 mm ID	N8122358
PFA/Sapphire Cassette	1.8 mm ID	N8122411
PFA/Platinum Torch Injector Assembly		
PFA/Platinum Twist Type	2.0 mm ID	N8122359
PFA/Platinum Cassette	2.0 mm ID	N8122412



Part No. **W1012406**



Part No. **N8122007**



Part No. **WE023951**

Injector Support Adapters

ICP-MS Model	Part No.
ELAN DRC/DRCplus/DRC II	
Ball Joint Twist Type Torch Injector Support Adapter	WE023951
ELAN DRC II	
Ball Joint Cassette Torch Injector Support Adapter	W1012406
ELAN 9000/6x00/5000/DRC-e	
Twist Type Torch Injector Support Adapter	N8122007
Twist Type Torch Injector Support Adapter (Assembled with O-Rings)	N8120116
ELAN 9000/6x00/5000/DRCs/NexION 300/350	
Injector Support Adapter O-Ring Kit	N8120100
Cassette Torch Injector Support Adapter For Non-Ball Joint Injectors	W1013266
Ball Joint Cassette Torch Injector Support Adapter	W1012406

Injector Support Adapter O-Rings

Position	Torch Type	Part No.
Internal	Twist/Cassette	09210011
External (Large)	Cassette	W1013545
External	Twist/Cassette	09210012
Internal	Cassette (Non-Ball Joint Only)	N8122010
Internal	Twist (Non-Ball Joint)	WE013060

Injector Adapter Kits and O-Rings

Description	Qty.	Part No.
For NexION 300/350 Injector Adapter Kit		
Kit Includes:		N8140360
Adapter	1	W1012406
O-Ring 3/8 in. ID x 1/16 in.	1	W1008572
O-Ring 1/4 in. ID x 1/16 in.	1	W1008573
Injector Spacer	1	W1012409
Injector Nut	1	W1012408
For NexION 300/350 Injector Adapter Kit 2		
Kit Includes:		N8140527
Injector O-Ring 0.239 ID x 0.070 WD	1	09210011
O-Ring 0.551 ID x 0.070 WD	1	N8122010
O-Ring 3/8 in. ID x 1/16 in.	1	W1008572
Adapter	1	W1012406
Injector Nut	1	W1012408
Injector Spacer	1	W1012409
O-Ring 5/8 in. ID x 1/8 in. WD	1	W1013545
For ELAN Cassette O-Ring Kit for Injector Support		
Kit Includes:		N8120100
Injector O-Ring 0.239 ID x 0.070 WD	2	09210011
Torch O-Ring 0.364 ID x 0.070 WD	2	09210012
O-Ring 0.551 ID x 0.070 WD	1	N8122010
O-Ring 5/8 in. ID x 1/8 in. WD	1	W1013545

Load Coils



LumiCoil™ Plasma Load Coil

ICP-MS Model	Part No.
NexION 1000/2000	N8150279



RF Load Coil Assemblies

ICP-MS Model	Part No.
ELAN 9000/6x00/DRCs/NexION 300/350	WE021816
ELAN 5000	WE012411