



## NexION 300/350Q/X/D Consumable Kits

Typical one year's supply of consumables.

Description	Qty.	Part No.
<b>NexION 300/350Q/X/D Consumables Kit</b> Kit Includes:		<b>N8140513</b>
Demountable Quartz Torch	2	<b>N8122006</b>
Load Coil	1	<b>WE021816</b>
Quartz Ball Joint Injector 2.0 mm ID	2	<b>WE023948</b>
Cassette Torch Injector Support Adapter for Ball Joint Injectors	1	<b>W1012406</b>
Glass Cyclonic Spray Chamber	1	<b>N8145014</b>
Nickel Skimmer Cone	2	<b>W1026356</b>
Nickel Sampler Cone	2	<b>W1033612</b>
Hyper Skimmer Cone	1	<b>W1033995</b>
Sampler Cone Gasket	3	<b>W1040148</b>
Hyper Skimmer O-Ring	2	<b>09902123</b>
Hyper Skimmer Screw	2	<b>09919737</b>
Santoprene® 1.30 mm Gray/Gray Pump Tubing	4	<b>N0777444</b>
Flared PVC 0.38 mm Green/Orange Pump Tubing	4	<b>N0777042</b>
Drain Tubing (ft.)	6	<b>02506495</b>
PTFE Sample Uptake Tubing (cm)	12	<b>WE024375</b>
Instrument Filter (left)	2	<b>W1036712</b>
Instrument Filter (back right)	2	<b>W1036713</b>
Spray Chamber Drain Fitting	2	<b>N8145015</b>
Nebulizer Fitting	2	<b>N8145016</b>
Male Barb CTFE Fitting	2	<b>N8145017</b>
Female Barb CTFE Fitting	2	<b>N8145018</b>

## NexION 300/350S Consumable Kits

Typical one year's supply of consumables.

Description	Qty.	Part No.
<b>NexION 300/350S Consumables Kit</b> Kit Includes:		<b>N8140514</b>
Demountable Quartz Torch	2	<b>N8122006</b>
Load Coil	1	<b>WE021816</b>
Quartz Ball Joint Injector 2.0 mm ID	2	<b>WE023948</b>
Cassette Torch Injector Support Adapter for Ball Joint Injectors	1	<b>W1012406</b>
Quartz Cyclonic Spray Chamber	1	<b>N8145013</b>
Platinum Skimmer Cone	1	<b>W1026907</b>
Platinum Sampler Cone	1	<b>W1033614</b>
Hyper Skimmer Cone	1	<b>W1033995</b>
Sampler Cone Gasket	3	<b>W1040148</b>
Hyper Skimmer O-Ring	2	<b>09902123</b>
Hyper Skimmer Screw	2	<b>09919737</b>
Santoprene® 1.30 mm Gray/Gray Pump Tubing	4	<b>N0777444</b>
Flared PVC 0.38 mm Green/Orange Pump Tubing	4	<b>N0777042</b>
Drain Tubing (ft.)	6	<b>02506495</b>
PTFE Sample Uptake Tubing (cm)	12	<b>WE024375</b>
Instrument Filter (left)	1	<b>W1036712</b>
Instrument Filter (back right)	1	<b>W1036713</b>
Spray Chamber Drain Fitting	2	<b>N8145015</b>
Nebulizer Fitting	2	<b>N8145016</b>
Male Barb CTFE Fitting	2	<b>N8145017</b>
Female Barb CTFE Fitting	2	<b>N8145018</b>