



TUBING

QLA Part #	Description	OEM Cross Ref#
PPTUBV-01	Peristaltic Pump Tubing, Santoprene, coded purple and white (SET/8)	17-6000
PPTUBV-02	Peristaltic Pump Tubing, Silicone, coded purple and white (SET/8)	17-6010
PPTUBD-01	Peristaltic Pump Tubing, Santoprene (SET/6)	5640-0337
TAATUB-HP	3 Collar Tygon Auto Analyzer Tubing, 2.06mm ID (SET/12)	5041-2166
TEFTUB-HP	FEP Teflon Tubing for HP, 1/32" ID x 1/16" OD, (500' Roll)	QLA
TEFTUB-8V	Color-coded Sample Teflon Tubing, .042" ID x .066" OD, (8 Lengths of 50')	17-6035
TUBING	Tygon Tubing, 1/16" x 1/8" (Sold by Foot)	QLA



GLASS ACCESSORIES

QLA Part #	Description	OEM Cross Ref#
------------	-------------	----------------

VOLUMETRIC FLASKS

QLA's round bottom Volumetric Flasks are Class A calibrated to contain and to deliver. We offer both 20°C and 37°C flasks.

VOLFLK-05	500mL Volumetric Flask, Class A/20°C, Round Bottom	QLA
VOLFLK-09	900mL Volumetric Flask, Class A/20°C, Round Bottom	QLA
VOLFLK-1K	1000mL Volumetric Flask, Class A/20°C, Round Bottom	QLA
VOLFLK-X1	1000mL Volumetric Flask, Class A/37°C, Round Bottom	QLA
VOLFLK-X5	500mL Volumetric Flask, Class A/37°C, Round Bottom	QLA
VOLFLK-X9	900mL Volumetric Flask, Class A/37°C, Round Bottom	QLA

(Photos on next page)