

SilTite® FingerTite ferrules continued

Description	Column ID	Ferrule ID	Pack size	Part no.
SilTite FingerTite Shimadzu 2010 injector starter kit	0.53 mm	0.7 mm	*	07362054
SilTite FingerTite Shimadzu 17A INJ / FID starter kit	0.53 mm	0.7 mm	*	07362853
SilTite FingerTite Shimadzu 17A injector starter kit	0.53 mm	0.7 mm	*	07362854
Replacement parts				
SilTite FingerTite ferrule 0.4 mm	0.1-0.25 mm	0.4 mm	10	073630
SilTite FingerTite ferrule 0.5 mm	0.32 mm	0.5 mm	10	073631
SilTite FingerTite ferrule 0.7 mm	0.53 mm	0.7 mm	10	073632
SilTite FingerTite ferrule blanking	–	–	2	073633
SilTite FingerTite female nut	–	–	5	073636

* Each starter kit includes all the parts necessary to convert one GC system (one injector and one detector) to the SilTite FingerTite system. In addition there are five SilTite FingerTite nuts, ten SilTite FingerTite ferrules, and a ferrule install tool which allows you to seat the ferrule in the correct position on the capillary column.

Ferrules

Column ID	Description	Pack size	Part no.
GC17A, GC2010 and GC2014 detector / injectors (not for MS interfaces or QP2010 injector)			
0.10-0.32 mm ID columns	100% graphite	10	0726080
0.45-0.53 mm ID columns	100% graphite	10	0726082
5 mm OD packed columns	100% graphite	10	0726001
0.10-0.25 mm ID columns	SilTite metal - initial installation	10*	073350
0.10-0.25 mm ID columns	SilTite ferrules	10	073227
0.32 mm ID columns	SilTite metal - initial installation	10*	073351
0.32 mm ID columns	SilTite ferrules	10	073228
0.45-0.53 mm ID columns	SilTite metal - initial installation	10*	073352
0.53 mm ID columns	SilTite ferrules	10	073229
n/a	SilTite metal nuts - slotted	5	073232
QP5000/5050 standard MS interface			
QP5000-I 0.10-0.25 mm ID columns	15% graphite / 85% Vespel® ferrules	10	0726563
QP5000-I 0.32 mm ID columns	15% graphite / 85% Vespel® ferrules	10	0726564
QP5000-II & QP5050 0.10-0.25 mm ID columns	15% graphite / 85% Vespel® ferrules	10	0726561
QP5000-II & QP5050 0.32 mm ID columns	15% graphite / 85% Vespel® ferrules	10	0726562
0.10-0.25 mm ID columns	SilTite metal - initial installation	10*	073204
0.10-0.25 mm ID columns	SilTite ferrules	10	073227
0.32 mm ID columns	SilTite metal - initial installation	10*	073205
0.32 mm ID columns	SilTite ferrules	10	073228
0.32 mm ID columns	SilTite ferrules	10	073229
n/a	SilTite metal nuts - QP5000/5050 standard MS interface	5	073233
QP5000/5050 wide bore MS interface, QP2010 injector and QP2010 standard MS interface			
0.10-0.25 mm ID columns	15% graphite / 85% Vespel® ferrules	10	072663
0.32 mm ID columns	15% graphite / 85% Vespel® ferrules	10	072654
0.45-0.53 mm ID columns	15% graphite / 85% Vespel® ferrules	10	072655
0.10-0.25 mm ID columns	SilTite metal - initial installation	10*	073200
0.10-0.25 mm ID columns	SilTite ferrules	10	073220
0.32 mm ID columns	SilTite metal - initial installation	10*	073201
0.32 mm ID columns	SilTite ferrules	10	073221
0.45-0.53 mm ID columns	SilTite metal - initial installation	10*	073202
0.45-0.53 mm ID columns	SilTite ferrules	10	073222
n/a	SilTite metal nuts	5	073224
Replacement SilTite® nuts			
2010 GC-MS system		5	073224
2010 GC-MS system with QP5000 series MS		5	073224
2010/2014 GC injectors and atmospheric detectors		5	073224
QP5000 jet separator MS interface		5	073224
QP5000 direct MS interface		5	073233
All injectors jet separator (starter kit), except 2010/2014		5	073232

* Includes ten ferrules, two SilTite nuts.