

## Ferrules continued

Instrument	Column ID	Ferrule ID	Pack size	Part no.
<b>15% graphite / 85% Vespel® ferrules</b>				
For TRACE 1300 injectors and detectors at atmospheric pressure e.g. FID	0.1-0.25 mm	0.4 mm	10	073109
	0.32 mm	0.5 mm	10	073111
	0.53 mm	0.8 mm	10	073113
	for 1/8" OD packed columns	1/8"	10	072669
	for 1/4" OD packed columns	1/4"	10	072667
For TRACE 1300 GC-MS interface connection	0.1-0.25 mm	0.4 mm	10	072663
	0.32 mm	0.5 mm	10	072654
	0.53 mm	0.8 mm	10	072655
<b>100% graphite ferrules</b>				
For TRACE 1300 injectors and detectors at atmospheric pressure e.g. FID (not for GC-MS)	0.1-0.32 mm	0.5 mm	10	072635
	0.45-0.53 mm	0.8 mm	10	072636
	for 1/8" OD packed columns	1/8"	10	072602
	for 1/4" OD packed columns	1/4"	10	072601
<b>Replacement SilTite® metal ferrules</b>				
For All GC-MS interface connections	0.1-0.25 mm	0.4 mm	10	073330
	0.32 mm	0.5 mm	10	073331
	0.53 mm	0.8 mm	10	073332
	1/32"	0.81 mm	10	073333
For TRACE 1300 connections	0.1-0.25 mm	0.4 mm	10	073220
	0.32 mm	0.5 mm	10	073221
	0.53 mm	0.8 mm	10	073222
	1/32"	0.81 mm	10	073219
<b>Replacement SilTite® nuts</b>				
SilTite metal nuts	–	–	5	073230
For TRACE 1300 GC-MS interface connection	–	–	5	073224
For TRACE 1300 split / splitless injector	–	–	5	073226
<b>Replacement SilTite® base seals</b>				
For TRACE 1300 split / splitless injector	–	–	2	073400
	–	–	10	073401

\* Includes ten ferrules, two SilTite nuts. # To be used in combination with brass nut (part no. 1034085).