

Ferrules

Instrument	Column ID	Ferrule ID	Pack size	Part no.
15% graphite / 85% Vespel® ferrules				
Injectors and detectors at atmospheric pressure e.g. FID	0.1-0.25 mm	0.4 mm	10	073109
	0.32 mm	0.5 mm	10	073111
	0.53 mm	0.8 mm	10	073113
	for 1/8" OD packed columns	1/8"	10	072669
	for 1/4" OD packed columns	1/4"	10	072667
GC-MS interface connection	0.1-0.25 mm	0.4 mm	10	072663
	0.32 mm	0.5 mm	10	072654
	0.53 mm	0.8 mm	10	072655
100% graphite ferrules				
Injectors and detectors at atmospheric pressure e.g. FID (not for GC-MS)	0.1-0.32 mm	0.5 mm	10	072635
	0.45-0.53 mm	0.8 mm	10	072636
	for 1/8" OD packed columns	1/8"	10	072602
	for 1/4" OD packed columns	1/4"	10	072601
SilTite® metal ferrules				
GC-MS interface connection (starter kit)	0.1-0.25 mm	0.4 mm	10*	073200
	0.32 mm	0.5 mm	10*	073201
	0.53 mm	0.8 mm	10*	073202
Split / splitless injectors (starter kit)	0.1-0.25 mm	0.4 mm	10#	073270
	0.32 mm	0.5 mm	10#	073271
	0.45-0.53 mm	0.8 mm	10#	073272
	1/32"	0.81 mm	10#	073273
Replacement SilTite® metal ferrules				
All connections	0.1-0.25 mm	0.4 mm	10	073220
	0.32 mm	0.5 mm	10	073221
	0.53 mm	0.8 mm	10	073222
	1/32"	0.81 mm	10	073219
Replacement SilTite® nuts				
GC-MS interface connection	-	-	5	073224
Split / splitless injector	-	-	5	073226
Replacement SilTite® base seals				
Split / splitless injector	-	-	2	073400
	-	-	10	073401

* Includes ten ferrules, two SilTite nuts. # Includes ten ferrules, two SilTite nuts and two SilTite inlet base seals.