

## SilTite® FingerTite ferrules



Description	Column ID	Ferrule ID	Pack size	Part no.
SilTite FingerTite PerkinElmer injector/GC-MS starter kit	0.1-0.25 mm	0.4 mm	*	073623
SilTite FingerTite PerkinElmer injector/FID starter kit	0.1-0.25 mm	0.4 mm	*	073622
SilTite FingerTite PerkinElmer injector/FID starter kit	0.53 mm	0.7 mm	*	07362253
SilTite FingerTite PerkinElmer injector starter kit	0.53 mm	0.7 mm	*	07362254
<b>Replacement parts</b>				
SilTite FingerTite ferrule 0.4 mm	0.1-0.25 mm	0.4 mm	10	073630
SilTite FingerTite ferrule 0.5 mm	0.32 mm	0.5 mm	10	073631
SilTite FingerTite ferrule 0.7 mm	0.53 mm	0.7 mm	10	073632
SilTite FingerTite blanking ferrule	–	–	2	073633
SilTite FingerTite female nut	–	–	5	073636

\* Each starter kit includes all the parts necessary to convert one GC system (one injector and one detector) to the SilTite FingerTite system. In addition there are five SilTite FingerTite nuts, ten SilTite FingerTite ferrules, and a ferrule install tool which allows you to seat the ferrule in the correct position on the capillary column.

## Ferrules

Instrument	Column ID	Size of nut	Ferrule ID	Pack size	Part no.
<b>15% graphite / 85% Vespel® ferrules</b>					
For injectors and detectors at atmospheric pressure e.g. FID	0.1-0.25 mm	1/16"	0.4 mm	10	072663
	0.1-0.25 mm	1/8"	0.4 mm	10	0726703
	0.32 mm	1/16"	0.5 mm	10	072654
	0.32 mm	1/8"	0.5 mm	10	0726702
	0.45-0.53 mm	1/16"	0.8 mm	10	072655
	0.45-0.53 mm	1/8"	0.8 mm	10	072671
	for 1/8" OD packed columns	1/8"	1/8"	10	072669
	for 1/4" OD packed columns	1/4"	1/4"	10	072667
<b>100% graphite ferrules</b>					
Injectors and detectors at atmospheric pressure e.g. FID (not for GC-MS)	0.1-0.32 mm	1/16"	0.5 mm	10	072627
	0.1-0.32 mm	1/8"	0.5 mm	10	072624
	0.45-0.53 mm	1/16"	0.8 mm	10	072626
	0.45-0.53 mm	1/8"	0.8 mm	10	0726280
	1/8" OD packed columns	1/8"	1/8"	10	072622
	1/4" OD packed columns	1/4"	1/4"	10	072621
<b>SilTite® metal ferrules</b>					
GC-MS interface connection (starter kit)	0.1-0.25 mm	–	0.4 mm	10*	073200
	0.32 mm	–	0.5 mm	10*	073201
	0.53 mm	–	0.8 mm	10*	073202
<b>Replacement SilTite® ferrules</b>					
GC-MS interface connection	0.1-0.25 mm	–	0.4 mm	10	073220
	0.32 mm	–	0.5 mm	10	073221
	0.53 mm	–	0.8 mm	10	073222
	1/32"	–	0.81 mm	10	073219
<b>Replacement SilTite® nuts</b>					
SilTite metal nuts	–	–	–	5	073224

\* Includes ten ferrules, two SilTite nuts.