

O-rings and sealing rings

Description	Usage	Pack size	Part no.
For Thermo Scientific TRACE 1300 Series GC			
Viton o-ring	Temperatures up to 300°C. Suitable for liners with OD of 6.3 mm	10	0726532
Graphite sealing ring	Temperatures up to 450°C. Suitable for all inlet liners above except 092004 and 09200401	10	0726005
Graphite sealing ring	Temperatures up to 450°C. Suitable for use with liners 092004 and 09200401	10	0726006

SilTite® FingerTite ferrules



Description	Column ID	Ferrule ID	Pack size	Part no.
For Thermo Scientific TRACE 1300 Series GC				
SilTite FingerTite INJ / FID starter kit	0.1-0.25 mm	0.4 mm	*	073610
SilTite FingerTite capillary / FID starter kit	0.1-0.25 mm	0.4 mm	*	073611
SilTite FingerTite INJ / MS starter kit	0.1-0.25 mm	0.4 mm	*	073612
Replacement parts				
SilTite FingerTite ferrule 0.4 mm	0.1-0.25 mm	0.4 mm	10	073630
SilTite FingerTite ferrule 0.5 mm	0.32 mm	0.5 mm	10	073631
SilTite FingerTite ferrule 0.7 mm	0.53 mm	0.7 mm	10	073632
SilTite FingerTite ferrule blanking	–	–	2	073633
SilTite FingerTite female nut	–	–	5	073636
SilTite FingerTite INJ base seal	0.1-0.25 mm	–	2	073640
SilTite FingerTite MS adapter	–	–	1	0736102
SilTite FingerTite injector	0.1-0.25 mm	–	1	0736104

* Each starter kit includes all the parts necessary to convert one GC system (one injector and one detector) to the SilTite FingerTite system. In addition there are five SilTite FingerTite nuts, ten SilTite FingerTite ferrules, and a ferrule install tool which allows you to seat the ferrule in the correct position on the capillary column.

Ferrules

Instrument	Column ID	Size of nut	Ferrule ID	Pack size	Part no.
SilTite® metal ferrules					
GC-MS interface connection (starter kit)	0.1-0.25 mm	–	0.4 mm	10*	073450
	0.32 mm	–	0.5 mm	10*	073451
	0.53 mm	–	0.8 mm	10*	073452
For TRACE 1300 split / splitless injectors (starter kit)	0.1-0.25 mm	–	0.4 mm	10	073270
	0.32 mm	–	0.5 mm	10	073271
	0.45-0.53 mm	–	0.8 mm	10	073272
	1/32"	–	0.81 mm	10	073273
Replacement SilTite® nuts					
SilTite metal nuts	–	–	–	5	073224

* Includes ten ferrules, two SilTite nuts.