


5 µL to 10 µL guided plunger syringes

| Specifications | | | |
|---|------------------------|--------------------------------------|-------------------------|
| Accuracy and reproducibility | ±1% (dispensed volume) | Scale length | 5 µL to 10 µL = 54.1 mm |
| Borosilicate glass barrel outer diameter (OD) | 5 µL to 10 µL = 6.5 mm | International standards traceability | ✓ |

Guided plunger syringes are the most rugged syringe available, making them perfect for industrial environments.



5 µL to 10 µL guided plunger syringes

| Part number | Part description and detail | Replacement needle | Replacement plunger |
|-------------|---|--------------------|---------------------|
| 001400 | 5 µL fixed needle guided plunger syringe with 5 cm 0.47 mm OD bevel tipped needle | - | - |
| 001450 | 5 µL removable needle guided plunger syringe with 5 cm 0.47 mm OD bevel tipped needle | 036110 | - |
| 001495 | 5 µL removable needle guided plunger syringe with 5 cm 0.47 mm OD bevel tipped needle (half scale) | 037110 | - |
| 002400 | 10 µL fixed needle guided plunger syringe with 5 cm 0.47 mm OD bevel tipped needle | - | - |
| 002450 | 10 µL removable needle guided plunger syringe with 5 cm 0.47 mm OD bevel tipped needle | 037110 | - |
| 002455 | 10 µL removable needle guided plunger syringe with GT plunger and 5 cm 0.47 mm OD bevel tipped needle | 037110 | 031805 |