


## ■ Glass Filter Holders; 30, 250 ml

For Particle Counting





### Glass Filter Holders

Two compact vacuum filter holders for easy particulate analysis. Both the top and bottom part of the filter holders are easily and securely fastened together using the

metal clamp. The centering rim on the filter support ensures correct positioning of the membrane filter. The glass frit filter support guarantees uniform distribution of retained particles on the filter surface.

### □ Ordering Information

Description		Order No.
Glass filter holder	30 ml	16306
Filter diameter	25 mm (or 24 mm)	
Filtration area	3 cm <sup>2</sup>	
Capacity	30 ml	
Outlet	12 mm outer diameter	
Glass filter holder	250 ml	16307
Filter diameter	47 mm (or 50 mm)	
Filtration area	12.5 cm <sup>2</sup>	
Capacity	250 ml	
Outlet	15 mm outer diameter	



### Adapter, 16836 | Adapter, 16837

For use of a glass filter holder, 16306 or 16307, on a Combisart® stainless steel manifold.

### □ Ordering Information

Description	Order No.
Adapter with 11 mm opening in stopper; for using filter holder 16306 on a Combisart® manifold	16836
Replacement stopper for 16836	00280
Adapter with 14 mm opening in stopper; for using filter holder 16307 on a Combisart® manifold	16837
Replacement stopper for 16837	00281