


## ■ Polycarbonate Filter Holders

For Particle Counting




### Polycarbonate Filter Holder, 250 ml

This reusable, practical filter holder made of autoclavable plastic is ideal for analytical testing outside the laboratory. For use with 47 mm membrane filters.

Outlet: TR 20×2 mm male thread

### □ Ordering Information



#### Description

#### Order No.

Polycarbonate filter holder without receiver flask

16511

Polycarbonate filter holder with receiver flask

16510

Hand vacuum pump with gauge and 60 cm PVC tubing

16672