

Test equipment monitoring	<ul style="list-style-type: none"> – Internal calibration weight – Automatic reproTEST for the weighing system – Automatic temperature calibration with optional temperature adjustment set, YTM03MA – Documentation in accordance with the applicable quality assurance requirements
Safety	Parameter settings password-protected against unauthorized access
Voltage frequency	48 to 60 Hz
Power consumption	Max. 700 VA
Operating temperature range	10°C to 30°C
Housing dimensions (W × D × H)	350 mm × 453 mm × 156 mm
Weight	Approx. 8.0 kg

Ordering Information

Accessories

6965542	Disposable sample pans, aluminum, round (80/box)
6906940	Glass fiber pads, for moisture analysis of liquid, pasty, greasy or fatty samples (80/box), hard quality
6906941	Glass fiber pads, for moisture analysis of liquid, pasty, greasy or fatty samples (200/box), soft quality
YCW4528-02	External calibration weight, 1 × 50 g, E2 with DAkkS certificate*
YDP01MA	Internal data printer, integrated into the device
YDP20-OCE	Data printer, external, verified for use in legal metrology; with date, time, statistics and transaction counter functions
6906918	Ink ribbon cassette, for YDP20-OCE and YDP01MA
6906937	Paper rolls for data printer, 5 rolls, 50 m each

* formerly known as DKD; DAkkS is the Germany's national accreditation body; DAkkS calibration certificates are internationally recognized