

HPLC Accessories

PEEK HPLC Fittings: PEEK Nuts & Ferrules

Upchurch PEEK Fingertight Connections

Upchurch PEEK Fingertight Union

Operate the PEEK zero dead volume unions to 5,000 psi (352 kg/cm²), the PEEK tees and crosses to 4,000 psi (281 kg/cm²). All items include Fingertight III nuts and ferrules.



Description	Cat. No.	Qty
Union, bore 0.010 in., thread: 10-32	57658	1 ea
Union, bore 0.020 in., thread: 10-32	57659	1 ea
Union, bore 0.020 in., thread: 1/4-28	57660-U	1 ea
Cross, bore 0.020 in., thread: 10-32	57663	1 ea
Tee, bore 0.020 in., thread: 10-32	57661	1 ea
Tee, bore 0.020 in., thread: 1/4-28	57662-U	1 ea

RheFlex® PEEK Fittings

For Rheodyne valves. Hold to 5000 psi/350kg/cm².



Long fitting set (57690-U)

Description	Cat. No.	Qty
Short fitting set, configured for 1/16 in. tubing	57691	5 ea
Long fitting set, configured for 1/16 in. tubing	57690-U	5 ea
Ferrule pack, configured for 1/16 in. tubing	57692	5 ea
Fitting set, configured for 1/8 in. tubing	57477	1 ea

Stainless Steel Fittings

Rheodyne® Fittings



from left to right: 58258, 58256, 58257

Description	Cat. No.	Qty
long nut, for 1/16 in. tubing	58257	10 ea
ferrule, for 1/16 in. tubing	58258	10 ea
nut, for 1/16 in. tubing	58256	10 ea
nut and ferrule, for 1/8 in. tubing	57479	1 ea

Each connection requires one ferrule and one male nut.

Valco® Fittings



From left to right: 22990-U, 22989

Description	Cat. No.	Qty
ferrule, configured for 1/16 in. tubing	22988	10 ea
ferrule, configured for 1/8 in. tubing	22989	10 ea
ferrule, for 1/4 in. tubing	58245-U	10 ea
nut, for for 1/16 in. tubing	22990-U	10 ea
nut, for for 1/8 in. tubing	22991	10 ea

Waters® fittings



58458

Description	Cat. No.	Qty
1/16 in.	58458	5 ea

SSI™ Fittings



From left to right: 58760-U, 58766

Description	Cat. No.	Qty
cap: 1/4-28, stainless steel	58766	1 ea
for 1/16 in. tubing	58764	1 ea
nut: 1/4-28	58760-U	10 ea
ferrule, for for 1/16 in. tubing	58762	10 ea