

## GC Accessories

Inlet Liners, Glass Wool, and Specialized Hand Tools: *Inlet Liners for Agilent (5890, 6890, and 7890)***Inlet Liner, Split/Splitless Type, Straight FocusLiner™ Design (wool packed)**

L × O.D. × I.D. .... 78.5 mm × 6.3 mm × 4.0 mm



Cat. No.	Qty
2879801-U	1 ea
2879805-U	5 ea
2879825-U	25 ea

**Inlet Liner, Split/Splitless Type, Straight Fast FocusLiner™ Design (wool packed)**

L × O.D. × I.D. .... 78.5 mm × 6.3 mm × 2.3 mm



Cat. No.	Qty
2879601-U	1 ea
2879605-U	5 ea
2879625-U	25 ea

**Inlet Liner, Split/Splitless Type, Straight Design (wool packed)**

L × O.D. × I.D. .... 78.5 mm × 6.3 mm × 4.0 mm



Cat. No.	Qty
2048601	1 ea
2048605	5 ea
2048625	25 ea

**Inlet Liner, Splitless Type, Single Taper Design (unpacked)**

L × O.D. × I.D. .... 78.5 mm × 6.5 mm × 4.0 mm



Cat. No.	Qty
2046601	1 ea
2046605	5 ea
2046625	25 ea

**Inlet Liner, Splitless Type, Dual-Taper Design (unpacked)**


L × O.D. × I.D. .... 78.5 mm × 6.5 mm × 4.0 mm



Cat. No.	Qty
2048501	1 ea
2048505	5 ea
2048525	25 ea

**Inlet Liner, Splitless Type, Recessed Gooseneck Design (unpacked)**


L × O.D. × I.D. .... 78.5 mm × 6.3 mm × 2.0 mm



Cat. No.	Qty
2879301-U	1 ea
2879305-U	5 ea

**Inlet Liner, Splitless Type, Straight Design (unpacked)**

L × O.D. × I.D. .... 78.5 mm × 6.3 mm × 4.0 mm



Cat. No.	Qty
2879401-U	1 ea
2879405-U	5 ea
2879425-U	25 ea

**Inlet Liner, Splitless Type, Straight Design (unpacked)**

L × O.D. × I.D. .... 78.5 mm × 6.5 mm × 2.0 mm



Cat. No.	Qty
2051301	1 ea
2051305	5 ea
2051325	25 ea

**Inlet Liner, Direct Type, Straight Design (unpacked)**

L × O.D. × I.D. .... 78.5 mm × 6.3 mm × 1.5 mm



Cat. No.	Qty
2051701	1 ea
2051705	5 ea
2051725	25 ea

**Inlet Liner, Direct (SPME) Type, Straight Design (unpacked)**

L × O.D. × I.D. .... 78.5 mm × 6.5 mm × 0.75 mm



Cat. No.	Qty
2637501	1 ea
2637505	5 ea
2637525	25 ea

**Inlet Liner, for Packed GC (unpacked)**

Deactivated (silanized) glass inlet liner for packed injection ports.

• Fits Agilent (5890, 6890, and 7890)

L × O.D. × I.D. .... 91.5 mm × 3.0 mm × 1.8 mm



Cat. No.	Qty
20508	5 ea
20511	25 ea