

GC Accessories

Inlet Liners, Glass Wool, and Specialized Hand Tools: *Inlet Liners for ATAS (Optic 2)**Inlet Liners for ATAS (Optic 2)*

Deactivated (silanized) glass inlet liners.

Inlet Liner, Split Type, Frit Design (unpacked)

L x O.D. x I.D. 80 mm x 5.0 mm x 3.0 mm




Cat. No.	Qty
2632605	5 ea
2632625	25 ea

Inlet Liners for Carlo Erba/Fisons (6000)

Deactivated (silanized) glass inlet liners.

Inlet Liner, Splitless Type, Straight Design (unpacked)

L x O.D. x I.D. 79.5 mm x 5.5 mm x 4.0 mm



Cat. No.	Qty
2632105	5 ea

Inlet Liner, Splitless Type, Straight Design (unpacked)

L x O.D. x I.D. 79.5 mm x 5.0 mm x 2.0 mm



Cat. No.	Qty
2632005	5 ea

Inlet Liners for Finnigan (9001GCQ)

Deactivated (silanized) glass inlet liners.

Inlet Liner, Split Type, Cup Design (wool packed)

L x O.D. x I.D. 78.5 mm x 6.3 mm x 4.0 mm



Cat. No.	Qty
2048201	1 ea
2048205	5 ea
2048225	25 ea

Inlet Liner, Split Type, Cup Design (packed with 10% OV-1 on Chromosorb W HP)

L x O.D. x I.D. 78.5 mm x 6.3 mm x 4.0 mm



Cat. No.	Qty
2055101	1 ea
2055105	5 ea
2055125	25 ea

Inlet Liner, Split Type, Cup Design (unpacked)

L x O.D. x I.D. 78.5 mm x 6.3 mm x 4.0 mm

Cat. No.	Qty
2051001	1 ea
2051005	5 ea
2051025	25 ea

Inlet Liner, Split/Splitless Type, Single Taper FocusLiner™ Design (wool packed)


L x O.D. x I.D. 78.5 mm x 6.3 mm x 4.0 mm



Cat. No.	Qty
2879901-U	1 ea
2879905-U	5 ea
2879925-U	25 ea

Inlet Liner, Split/Splitless Type, Single Taper Fast FocusLiner™ Design (wool packed)

L x O.D. x I.D. 78.5 mm x 6.3 mm x 2.3 mm



Cat. No.	Qty
2879501-U	1 ea
2879505-U	5 ea
2879525-U	25 ea

Inlet Liner, Split/Splitless Type, Single Taper Design (wool packed)


L x O.D. x I.D. 78.5 mm x 6.5 mm x 4.0 mm



Cat. No.	Qty
2047801	1 ea
2047805	5 ea
2047825	25 ea

Inlet Liner, Split/Splitless Type, Recessed Gooseneck Design (wool packed)

L x O.D. x I.D. 78.5 mm x 6.3 mm x 4.0 mm



Cat. No.	Qty
2879701-U	1 ea
2879705-U	5 ea
2879725-U	25 ea

Inlet Liner, Split/Splitless Type, Straight FocusLiner™ Design (wool packed)

L x O.D. x I.D. 78.5 mm x 6.3 mm x 4.0 mm



Cat. No.	Qty
2879801-U	1 ea
2879805-U	5 ea
2879825-U	25 ea