

GC Accessories

Inlet Liners, Glass Wool, and Specialized Hand Tools: *Inlet Liners for Varian® (CP-1177 Injector)***Inlet Liner, Split/Splitless Type, Single Taper FocusLiner™ Design (wool packed)**

L × O.D. × I.D. 78.5 mm × 6.3 mm × 4.0 mm



Cat. No.	Qty
2879901-U	1 ea
2879905-U	5 ea
2879925-U	25 ea

Inlet Liner, Split/Splitless Type, Single Taper Fast FocusLiner™ Design (wool packed)

L × O.D. × I.D. 78.5 mm × 6.3 mm × 2.3 mm



Cat. No.	Qty
2879501-U	1 ea
2879505-U	5 ea
2879525-U	25 ea

Inlet Liner, Split/Splitless Type, Single Taper Design (wool packed)


L × O.D. × I.D. 78.5 mm × 6.5 mm × 4.0 mm



Cat. No.	Qty
2047801	1 ea
2047805	5 ea
2047825	25 ea

Inlet Liner, Split/Splitless Type, Recessed Gooseneck Design (wool packed)

L × O.D. × I.D. 78.5 mm × 6.3 mm × 4.0 mm



Cat. No.	Qty
2879701-U	1 ea
2879705-U	5 ea
2879725-U	25 ea

Inlet Liner, Split/Splitless Type, Straight FocusLiner™ Design (wool packed)

L × O.D. × I.D. 78.5 mm × 6.3 mm × 4.0 mm



Cat. No.	Qty
2879801-U	1 ea
2879805-U	5 ea
2879825-U	25 ea

Inlet Liner, Split/Splitless Type, Straight Fast FocusLiner™ Design (wool packed)

L × O.D. × I.D. 78.5 mm × 6.3 mm × 2.3 mm



Cat. No.	Qty
2879601-U	1 ea
2879605-U	5 ea
2879625-U	25 ea

Inlet Liner, Split/Splitless Type, Straight Design (wool packed)

L × O.D. × I.D. 78.5 mm × 6.3 mm × 4.0 mm



Cat. No.	Qty
2048601	1 ea
2048605	5 ea
2048625	25 ea

Inlet Liner, Splitless Type, Single Taper Design (unpacked)

L × O.D. × I.D. 78.5 mm × 6.5 mm × 4.0 mm



Cat. No.	Qty
2046601	1 ea
2046605	5 ea
2046625	25 ea

Inlet Liner, Splitless Type, Dual-Taper Design (unpacked)


L × O.D. × I.D. 78.5 mm × 6.5 mm × 4.0 mm



Cat. No.	Qty
2048501	1 ea
2048505	5 ea
2048525	25 ea

Inlet Liner, Splitless Type, Recessed Gooseneck Design (unpacked)


L × O.D. × I.D. 78.5 mm × 6.3 mm × 2.0 mm



Cat. No.	Qty
2879301-U	1 ea
2879305-U	5 ea

Inlet Liner, Splitless Type, Straight Design (unpacked)

L × O.D. × I.D. 78.5 mm × 6.3 mm × 4.0 mm



Cat. No.	Qty
2879401-U	1 ea
2879405-U	5 ea
2879425-U	25 ea