


Colorimeter Cells

Colorimeter cells are fully fused glass cells for use with colorimeters. Each cell has two polished windows and an open top with no cover.



Catalog Number	Path Length mm	Exterior, mm			Interior, mm		Nominal Vol. ml	Polished Windows	Price per cell
		Width	Length	Height	Width	Length			
<i>For use with: Hunter Colorimeter, ACS, Data Color and others</i>									
93-G-2	2	55	2	56	51	2	25	2	\$ 215.50
93-G-5	5	55	9	56	51	5	12.5	2	173.20
93-G-10	10	55	14	56	51	10	25	2	156.60
93-G-20	20	55	24	56	51	20	50	2	193.40
93-G-30	30	55	34	56	51	30	75	2	269.60
93-G-50	50	55	54	56	51	50	125	2	313.20
<i>For use with: CP method and general colorimetric work</i>									
97-G-40H	40	20	44	45	16	40	25	2	\$ 85.10
<i>For use with: Absorptimeters</i>									
96-G-2.5	2.5	28	8.5	40	24	2.5	1.8	2	\$ 45.20
96-G-5	5	28	11	40	24	5	3.6	2	45.20
96-G-10	10	28	16	40	24	10	7.2	2	45.20
96-G-20	20	28	26	40	24	20	14	2	64.10
96-G-40	40	28	46	40	24	40	28	2	85.20
96-G-100	100	28	106	40	24	100	70	2	135.60