

General laboratory accessories

In addition to the filtration consumable range, we provide a comprehensive range of accessories for routine work in your laboratory. The table below shows a selection of the products we offer.




1PS phase separator



Grade 105 lens cleaning tissue



Benchkote™ protection paper



pH paper



Vacu-Guard Pump protection filter

Description	Product name	Dimension	Code no.	Quantity
Phase separation paper <ul style="list-style-type: none"> • Separatory funnel replacement: Automatic cut-off • Ease of use: No special training required 	1PS Phase separator paper	Diam. 125 mm	2200-125	100/pack
		Diam. 150 mm	2200-150	100/pack
Optical lens cleaning tissue <ul style="list-style-type: none"> • Soft tissue for removing surface moisture and grease from lenses and other optical surfaces 	Grade 105	100 × 150 mm	2105-841	25 wallets of 25 sheets
		200 × 300 mm	2105-862	100/pack
Benchkote bench protection papers <ul style="list-style-type: none"> • High-quality, smooth, absorbent Whatman paper • Quickly absorbs liquid spills and protect the working surface • Benchkote Plus is thicker and more absorbent 	Benchkote	460 × 570 mm	2300-916	50/pack
		460 mm × 50 m	2300-731	1/pack
	Benchkote Plus	500 × 600 mm	2301-6150	50/pack
		600 mm × 50 m	2301-6160	1/pack
pH Indicator Paper <ul style="list-style-type: none"> • Range of pH indicator and test papers for the rapid determination of pH values in many applications 	Color Bonded, 0.0 to 14.0 range	6 × 80 mm	2613-991	100 strips, 1/pack
	Standard Full Range, Reel, 1.0 to 14.0 range	7 mm × 5 m	2600-100A	1/pack
	Standard Narrow Range, Reel, 4.0 to 7.0 range	7 mm × 5 m	2600-102A	1/pack
Pump protection filters <ul style="list-style-type: none"> • Protects vacuum pump systems from aqueous aerosols. Hydrophobic PTFE membranes retain 99,99% of airborne particles > 0.1 µm 	Vacu-Guard	50 mm	6722-5000	10/pack