

Ordering Information

Oasis HLB Sample Extraction Products

| Description | Format | Particle Size | Qty. | P/N |
|---|---------------|---------------|----------|-------------|
| Oasis HLB Cartridge | 1 cc/10 mg | 30 µm | 100/box | 186000383 |
| Oasis HLB Cartridge | 1 cc/30 mg | 30 µm | 100/box | WAT094225 |
| Oasis HLB Cartridge | 1 cc/30 mg | 30 µm | 1000/box | 186003908 |
| Oasis HLB Flangeless Cartridge | 1 cc/30 mg | 30 µm | 100/box | 186001879 |
| Oasis HLB Cartridge with Gilson ASPEC Adapter | 1 cc/10 mg | 30 µm | 500/box | 186000988 |
| Oasis HLB Cartridge with Gilson ASPEC Adapter | 1 cc/30 mg | 30 µm | 500/box | WAT058882 |
| Oasis HLB Cartridge | 3 cc/60 mg | 30 µm | 100/box | WAT094226 |
| Oasis HLB Cartridge | 3 cc/60 mg | 30 µm | 1000/box | 186007646 |
| Oasis HLB Flangeless Cartridge | 3 cc/60 mg | 30 µm | 100/box | 186001880 |
| Oasis HLB Cartridge with Gilson ASPEC Adapter | 3 cc/60 mg | 30 µm | 500/box | WAT058883 |
| Oasis HLB Cartridge | 6 cc/200 mg | 30 µm | 30/box | WAT106202 |
| Oasis HLB Cartridge | 3 cc/400 mg | 60 µm | 100/box | 186003849 |
| Oasis HLB Cartridge | 3 cc/540 mg | 60 µm | 100/box | 186004134 |
| Oasis PRiME HLB Cartridge | 3 cc/150 mg | — | 100/box | 186008717 |
| Oasis PRiME HLB Cartridge | 3 cc/150 mg | — | 30/box | 186008718 |
| Oasis HLB Flangeless Cartridge | 3 cc/540 mg | 60 µm | 100/box | 186003852 |
| Oasis HLB Cartridge | 6 cc/150 mg | 30 µm | 30/box | 186003365 |
| Oasis HLB Cartridge | 6 cc/150 mg | 60 µm | 30/box | 186003379 |
| Oasis HLB Cartridge | 6 cc/500 mg | 60 µm | 30/box | 186000115 |
| Oasis HLB Cartridge | 12 cc/500 mg | 60 µm | 20/box | 186000116 |
| Oasis HLB Cartridge | 20 cc/1 g | 60 µm | 20/box | 186000117 |
| Oasis HLB Cartridge | 35 cc/6 g | 60 µm | 10/box | 186000118 |
| Oasis HLB Plus Short Cartridge | 225 mg | 60 µm | 50/box | 186000132 |
| Oasis HLB Plus Light Cartridge | 30 mg | 30 µm | 50/box | 186005125 |
| Oasis HLB Vac RC Cartridge | 20 cc/30 mg | 30 µm | 50/box | 186000382 |
| Oasis HLB Vac RC Cartridge | 20 cc/60 mg | 30 µm | 50/box | 186000381 |
| Oasis HLB Glass Cartridge | 5 cc/200 mg | 60 µm | 30/box | 186000683 |
| Oasis HLB µElution Plate | 2 mg/96-well | 30 µm | 1/pk | 186001828BA |
| Oasis HLB Plate | 5 mg/96-well | 30 µm | 1/pk | 186000309 |
| Oasis HLB Plate | 10 mg/96-well | 30 µm | 1/pk | 186000128 |
| Oasis HLB Plate | 30 mg/96-well | 30 µm | 1/pk | WAT058951 |
| Oasis HLB Plate | 60 mg/96-well | 60 µm | 1/pk | 186000679 |

OASIS MCX FOR BASIC COMPOUNDS

Obtain selective retention of basic drugs with cation-exchange groups on the sorbent surface. The Oasis MCX (Mixed-mode Cation eXchange) Sorbent has a tightly controlled ion-exchange capacity (1 meq/g). There are no silanol groups to complicate the retention mode or method development. This novel, water-wettable, polymeric sorbent is stable from pH 0–14, making method development simple and fast.

Ordering Information

Oasis MCX Sample Extraction Products (Cation Exchange)

| Description | Format | Particle Size | Qty. | P/N |
|--------------------------------|---------------|---------------|---------|-------------|
| Oasis MCX Cartridge | 1 cc/10 mg | 30 µm | 100/box | 186004648 |
| Oasis MCX Cartridge | 1 cc/30 mg | 30 µm | 100/box | 186000252 |
| Oasis MCX Flangeless Cartridge | 1 cc/30 mg | 30 µm | 100/box | 186001881 |
| Oasis MCX Cartridge | 1 cc/30 mg | 60 µm | 100/box | 186000782 |
| Oasis MCX Cartridge | 3 cc/60 mg | 30 µm | 100/box | 186000254 |
| Oasis MCX Flangeless Cartridge | 3 cc/60 mg | 30 µm | 100/box | 186001882 |
| Oasis MCX Cartridge | 3 cc/60 mg | 60 µm | 100/box | 186000253 |
| Oasis MCX Cartridge | 6 cc/150 mg | 30 µm | 30/box | 186000256 |
| Oasis MCX Cartridge | 6 cc/150 mg | 60 µm | 30/box | 186000255 |
| Oasis MCX Cartridge | 6 cc/500 mg | 60 µm | 30/box | 186000776 |
| Oasis MCX Cartridge | 20 cc/1 g | 60 µm | 20/box | 186000777 |
| Oasis MCX Cartridge | 35 cc/6 g | 60 µm | 10/box | 186000778 |
| Oasis MCX Plus Short Cartridge | 225 mg | 60 µm | 50/box | 186003516 |
| Oasis MCX Vac RC Cartridge | 20 cc/60 mg | 30 µm | 50/box | 186000261 |
| Oasis MCX Vac RC Cartridge | 20 cc/60 mg | 60 µm | 50/box | 186000380 |
| Oasis MCX µElution Plate | 2 mg/96-well | 30 µm | 1/pk | 186001830BA |
| Oasis MCX Plate | 10 mg/96-well | 30 µm | 1/pk | 186000259 |
| Oasis MCX Plate | 30 mg/96-well | 30 µm | 1/pk | 186000248 |
| Oasis MCX Plate | 30 mg/96-well | 60 µm | 1/pk | 186000250 |
| Oasis MCX Plate | 60 mg/96-well | 60 µm | 1/pk | 186000678 |

DID YOU KNOW...

Oasis Cartridges and Plates are available in two particle sizes (30 µm or 60 µm).

This allows you to select the appropriate product based on the viscosity and turbidity of your sample. For extraction of most plasma, serum, and human urine, choose the 30 µm sorbent. For more viscous samples such as animal urine, excellent flow can be achieved using the 60 µm sorbent in either cartridges or plates.