








AUTOSAMPLER VIALS FOR WATERS SYSTEMS

Vials for Waters 717 Autosampler

	4 mL Screw Neck	Amber Screw Neck	Total Recovery	PP Screw Neck Vial	PP Conical	Glass Shell Vial	Amber Glass Shell Vial
15 × 45 mm							
48-position Carousel	33	34	35	36	37	38	39

Combination Packs							
Vial, Cap, and LectraBond PTFE/Silicone Septum	186000838C	186001133C	186002629C	—	—	—	—
Combination Deactivated	186000838DV	186001133DV	—	—	—	—	—
Vial, Cap, and LectraBond Pre-slit PTFE/ Silicone Septum	186000839C	186001134C	186002630C	—	—	—	—
Combination Deactivated	186000839DV	186001134DV	—	—	—	—	—
Vial and PE Snap Cap	—	—	—	—	186004031	WAT025051	WAT025050

Components							
Vials Only	186000840(DV)	186001135(DV)	186002520	186000999 ¹	—	—	—
Max Volume Injection/Max Residual Volume	2400 µL/1600 µL	2400 µL/1600 µL	3000 µL/40 µL	2000 µL/1000 µL	2950 µL/50 µL	2400 µL/1600 µL	2400 µL/1600 µL
Cap LectraBond PTFE/Silicone 100/pk	186000841	186000841	186000841	—	—	—	—
Screw Cap with Bonded PTFE/Silicone Septum, 1000/pk	—	—	—	186000965	—	—	—
Cap LectraBond Pre-slit PTFE/Silicone, 100/pk	186000842	186000842	186000842	—	—	—	—
Black Phenol Cap, 144/pk	WAT072711	WAT072711	WAT072711	—	—	—	—
PTFE Septum, 1440/pk	WAT073005	WAT073005	WAT073005	—	—	—	—
PTFE Septum, 144/pk	WAT072714	WAT072714	WAT072714	—	—	—	—
Self Sealing Septum, 144/pk	WAT022861	WAT022861	WAT022861	—	—	—	—
250 µL Glass Insert ²	WAT072704(DV)	WAT072704(DV)	—	—	—	—	—
Max Volume Injection/Max Residual Volume	244 µL/6 µL	244 µL/6 µL	—	—	—	—	—
250 µL Glass Insert, 144/pk ²	WAT015199(DV)	WAT015199(DV)	—	—	—	—	—
Max Volume Injection/Max Residual Volume	230 µL/20 µL	230 µL/20 µL	—	—	—	—	—
250 µL Plastic Conical Insert (PMP), 144/pk ²	WAT072030	WAT072030	—	—	—	—	—
Max Volume Injection/Max Residual Volume	230 µL/20 µL	230 µL/20 µL	—	—	—	—	—
Springs for LVI, 100/pk	WAT072708	WAT072708	—	—	—	—	—

Storage Cap							
Solid Black Cap with Silicone/ PTFE Liner for Sample Storage	186007224	186007224	186007224	—	—	—	—

¹Item contains 1000 vials.

²Inserts requires springs, p/n: WAT072708.

When (DV) appears beside the part number, a deactivated version of this product can be ordered by adding DV to the right of the part number.