









Complete Listing of 12 × 32 mm Vials and Accessories

|  | Clear   | Amber   | Max Recovery  | 300 µL PP   | 750 µL PP  | 10 mm Cap Clear   | Total Recovery  | Amber Max   |
|--|---|---|---|---|--|---|---|---|
| 12 × 32 mm   |  |  |  |  |  |  |  |  |
| <b>Vial Number</b>   | 17  | 18  | 19  | 20  | 21   | 22  | 23  | 24  |
| <b>Compatible Systems</b>  |   |   |   |   |  |   |   |   |
| Alliance 2690/2695   | ▪   | ▪   | —   | ▪   | ▪  | ▪   | ▪   | —   |
| Alliance 2790/2795   | ▪   | ▪   | ▪   | ▪   | ▪  | —   | —   | ▪   |
| ACQUITY  | ▪   | ▪   | ▪   | ▪   | ▪  | —   | ▪   | ▪   |
| <b>Inserts</b>   |   |   |   |   |  |   |   |   |
| 300 µL with Poly Spring  | WAT094170(DV)   | WAT094170(DV)   | —   | —   | —  | WAT094170(DV)   | —   | —   |
| Max Volume Injection/<br>Max Residual Volume                         | 230 µL/20 µL  | 230 µL/20 µL  | —   | —   | —  | 230 µL/20 µL  | —   | —   |
| 150 µL with Poly Spring  | WAT094171(DV)   | WAT094171(DV)   | —   | —   | —  | WAT094171(DV)   | —   | —   |
| Max Volume Injection/<br>Max Residual Volume                         | 144 µL/6 µL   | 144 µL/6 µL   | —   | —   | —  | 144 µL/6 µL   | —   | —   |
| <b>Black Screw Cap for TruView Vials</b>                             |   |   |   |   |  |   |   |   |
| PTFE/Silicone Septum   | 186005826   | 186005826   | 186005826   | —   | —  | —   | 186005826   | 186005826   |
| Pre-slit PTFE/Silicone Septum  | 186005827   | 186005827   | 186005827   | —   | —  | —   | 186005827   | 186005827   |
| <b>Light Blue Screw Cap for LCMS Certified Vials</b>                 |   |   |   |   |  |   |   |   |
| PTFE/Silicone Septum   | 186005828   | 186005828   | 186005828   | —   | —  | —   | 186005828   | 186005828   |
| Pre-slit PTFE/Silicone Septum  | 186005829   | 186005829   | 186005829   | —   | —  | —   | 186005829   | 186005829   |
| <b>Screw Cap and Septum-Silicone/PTFE</b>                            |   |   |   |   |  |   |   |   |
| PE Septumless Cap  | 186004169   | 186004169   | 186004169   | 186004169   | 186004169  | —   | 186004169   | 186004169   |
| Blue LectraBond  | 186000274   | 186000274   | 186000274   | 186000274   | 186000274  | —   | 186000274   | 186000274   |
| Red LectraBond   | 186002129   | 186002129   | 186002129   | 186002129   | 186002129  | —   | 186002129   | 186002129   |
| Green LectraBond   | 186002130   | 186002130   | 186002130   | 186002130   | 186002130  | —   | 186002130   | 186002130   |
| White LectraBond   | 186002456   | 186002456   | 186002456   | 186002456   | 186002456  | —   | 186002456   | 186002456   |
| Black Cap with PTFE Septum, 100/pk                                   | 186007198   | 186007198   | 186007198   | 186007198   | 186007198  | —   | 186007198   | 186007198   |
| <b>Screw Cap and Pre-slit Septum-Silicone/PTFE</b>                   |   |   |   |   |  |   |   |   |
| Blue LectraBond  | 186000305   | 186000305   | 186000305   | 186000305   | 186000305  | —   | 186000305   | 186000305   |
| Red LectraBond   | 186002128   | 186002128   | 186002128   | 186002128   | 186002128  | —   | 186002128   | 186002128   |
| Green LectraBond   | 186002127   | 186002127   | 186002127   | 186002127   | 186002127  | —   | 186002127   | 186002127   |
| White LectraBond   | 186002457   | 186002457   | 186002457   | 186002457   | 186002457  | —   | 186002457   | 186002457   |
| <b>For Dissolution System</b>  |   |   |   |   |  |   |   |   |
| Pre-assembled Vial, Cap, and<br>Pre-slit Septum                      | 186000989(DV)   | 186003455   | —   | —   | —  | —   | —   | —   |
| <b>Storage Cap</b>   |   |   |   |   |  |   |   |   |
| Black Solid 9 mm Cap with Silicone/<br>PTFE Liner for Sample Storage | 186007187   | 186007187   | 186007187   | 186007187   | 186007187  | —   | 186007187   | 186007187   |

When (DV) appears beside a number, a deactivated version of the part can be ordered by adding a DV to the right of the part number.