











## Screw Cap Vials for Compatible Systems

	Clear	Amber	Amber Max Recovery	Clear Glass Max Recovery	Qsert Vial	Amber Qsert	PP 300 µL	PP 750 µL	10 mm Cap	PP 250 µL 8 mm Cap
12 × 32 mm										
<b>Vial Number</b>	<b>50</b>	<b>51</b>	<b>52</b>	<b>53</b>	<b>54</b>	<b>55</b>	<b>56</b>	<b>57</b>	<b>58</b>	<b>59</b>
<b>Compatible Systems</b>										
Agilent Technologies	▪	▪	▪	▪	▪	▪	▪	▪	—	—
Alcott, Antek, CTC, Spark Thermal Separations	—	—	—	—	—	—	—	—	▪	▪
Beckman, Dynatech, Finnigan, Fisons, Gilson	▪	▪	▪	▪	▪	▪	▪	▪	—	—
Hitachi, LDC, Perkin-Elmer, Shimadzu, Spectra-Physics, Thermo, Varian	▪	▪	▪	▪	▪	▪	▪	▪	▪	▪
<b>LCMS Certified Combination Packs</b>										
Vial, Cap, and Silicone/PTFE Septum	600000751CV	600000752CV	600000754CV	600000749CV	—	—	—	—	—	—
Vial, Cap, and Pre-slit Silicone/PTFE Septum	600000668CV	600000669CV	600000755CV	600000670CV	—	—	—	—	—	—
<b>LC/GC Certified Combination Packs</b>										
Bonded Silicone/PTFE Septum	186000272C	186000846C	186003885C	186000326C	186001126C	186001130C	—	—	WAT270946C <sup>1</sup>	—
Combination Deactivated <sup>2</sup>	186000272DV	186000846DV	—	186000326DV	186001126DV	186001130DV	—	—	WAT270946DV	—
Bonded Pre-slit Silicone/PTFE Septum	186000307C	186000847C	186003886C	186000327C	186001128C	186001131C	—	—	—	—
Combination Deactivated <sup>2</sup>	186000307DV	186000847DV	—	186000327DV	186001128DV	186001131DV	—	—	—	—
<b>Combination Packs</b>										
Bonded Silicone/PTFE Septum	—	—	—	—	—	—	186002640	186005220	—	—
Bonded Pre-slit Silicone/PTFE Septum	—	—	—	—	—	—	186002639	186005221	—	—
<b>Vials Only</b>										
Vials Only	186000273	186000848	—	186002802	186002804	186002803	186002626	186005219	WAT063300	WAT094172
Deactivated Vials Only	186000273DV	186000848DV	—	—	—	—	—	—	WAT063300DV	—
<b>Inserts</b>										
300 µL with Poly Spring	WAT094170	WAT094170	—	—	—	—	—	—	WAT094170	—
300 µL with Poly Spring Deactivated	WAT094170DV	WAT094170DV	—	—	—	—	—	—	WAT094170DV	—
150 µL with Poly Spring	WAT094171	WAT094171	—	—	—	—	—	—	WAT094171	—
150 µL with Poly Spring Deactivated	WAT094171DV	WAT094171DV	—	—	—	—	—	—	WAT094171DV	—

<sup>1</sup> Septum not bonded.

<sup>2</sup> Not certified.

All items come in quantities of 100 unless otherwise noted.